



RISK	High Risk / Elevated Risk
WATERLOW	Up to 20
WARRANTY	3 years
REPLACEMENTS	All parts
WEIGHT LIMIT	22 stone / 140 kg
HEIGHT	3" / 7.8 cm
COVERS	Two way stretch vapour permeable
PLATILON LINER	Standard
CUSHION SIZE RANGE	23 sizes from 14"x13" to 18"x16" / 35.6 x 33.1 cm to 45.7 x 40.6 cm
PRODUCT CODE	FI/Size



FOR HIGH RISK / ELEVATED RISK*

The Invacare Flo-tech Image cushion uses a combination of a slim, shaped foam cushion base and a larger gel sac increasing even weight distribution and immersion, thus delivering a high degree of protection to users who are at High Risk / Elevated Risk of developing pressure ulcers. This slim cushion format is ideal for both chair and wheelchair use.

FEATURES

Slimline contoured foam cushion

Gel sac offers protection to the ischial and sacral area

Moulded foam cushion

Pommel at front of the cushion

Contoured seating surface

Castellations within seating area allow slight movement

High quality resilient foam, manufactured using a unique formula

Integral sag compensator feature is available

All covers are two-way stretch, water resistant and vapour permeable

All covers have an anti-slip base

Handle on the base of the cushion

BENEFITS



Some users prefer a slimline cushion, finding it more compatible with their daily living activities

The coverage afforded by this large gel sac increases the support surface area, and reduces the occurrence of peak pressures

Moulding allows Flo-tech cushions to provide a 3D form, which cannot be emulated in cut foam

Separates the thighs, encourages good posture and allows air to circulate. Also helps to prevent the knees from rubbing together

Offers improved stability and postural support. Also delivers greater degree of weight distribution, compared to that of a flat surface

Reduces shear and friction forces

Offers high levels of pressure reduction, and product longevity

Designed to accommodate the sag in wheelchair seats

Allows patient movement – reduces potential for shear and friction and ensures good air circulation, reducing potential for skin maceration

Ensures safe use of the cushion and helps prevent slippage

Assists in the easy transfer of the cushion

*Essential nursing care is pivotal in pressure ulcer prevention. This cushion will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. The Waterlow score is a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (nursing judgement). Informal assessment is considered to be of greater importance and clinical value.