



Yes, you can.[®]

➤ *Invacare*[®] **Action⁴ Heavy Duty** **NEW**

Strength durability & easier handling

The *Invacare*[®] **Action⁴ Heavy Duty** manual chair has been developed to offer a durable and reliable wheelchair with a maximum user weight of 150 kg. It also offers compatibility with the *Invacare*[®] **Action³** range to make refurbishment and upgrades easy.



Easy rolling and manoeuvrability

The *Invacare® Action⁴ HD* has a reinforced chassis which provides a stable yet lightweight base that can be manoeuvred easily with minimum effort. It has a choice of seat widths from 50.5 cm to 60.5 cm and offers many adjustable features to optimize posture and comfort.




Compatibility

There are a number of different options available to customise this chair to client requirements, which includes a tension adjustable backrest, a range of multi-adjustable legrests, armrests, hub brakes, anti tip device, and optional rear wheel positions. All parts fitted to the frame of the *Action⁴ HD* are interchangeable with the *Invacare® Action³* to help keep down maintenance costs and to provide easier refurbishment.



**Adjustable Backrest
Upholstery**



Armrests



Hub brakes

Technical data

38, 40.5, 43, 45.5, 48, 50.5, 55.5, 60.5 cm	40-45 cm	38.5, 41, 43.5, 46, 48.5, 51 cm	38 - 46 cm 39 - 51 cm*	20, 23, 26, 19 - 29.5 cm	35 - 48 cm	0° - 6°	7° 0° - 30°	
seat width + 21 cm	90 - 103 cm	105 - 111 cm	30.5 cm	15.3 kg (Seat width 44 cm)	150 kg	Transport: Weight 9 kg	Crash tested ISO 7176-19**	

* Reclining backrest

** Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

Frame colour



Sapphire Blue



Anthracite Black



Pearl Grey



Black nylon

Please notice colors may vary slightly from those displayed above.

Upholstery

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited · South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016
E-mail: uk@invacare.com · www.invacare.co.uk

