



Pit Beirer
World Champion Motocross Runner Up

Neon



QUICKIE

INNOVATE.
BUILD FOR THE FUTURE

QUICKIE

Neon

The new Neon from Quickie re-defines Mini-Max principle. Minimum space requirement offering maximum performance.



Frame colour: bright silver
Sling: black

Small, smaller, Neon – a miracle of space



A clever frame



Automatically folding footrest

Beautifully slim

With its innovative cross-brace and the newly-designed, fold-down backrest, the Neon achieves folded dimensions unheard of up to now. With its "Dream dimensions" of (W x H x L) 25 cm x 49 cm x 78 cm, the new Neon is the smallest wheelchair in its class, and will fit in practically anywhere. Whether it's in the car or when travelling, at work or for leisure, the Neon is the ideal companion for the active person. Not only can the Neon make itself extremely small, it is also an absolute lightweight, with a weight of just 9.7 kg in its normal configuration.

The smallest fold-up!





Neon

Swing Away



Great looks, exceptional functionality

The new Neon is not only a trail-blazer in size and weight, but sets new standards in form and function. Its open, curved lines give the Neon an easy, elegant overall appearance. With Swing Away, the ability to swing the footrests inwards and outwards it offers maximum flexibility. The overall weight of only 10.5 kg means that the Neon Swing Away version is also extremely light.

QUICKIE

The Mini-Max-Principle



1



2



3



4

The new Neon is available in various widths, with continuous or two-piece footplates through to the newly-developed footrests which can swing inwards and outwards, and an automatically folding footrest.

The new Neon is the perfect chair for all requirements and is totally convincing in every aspect. That's what we mean by the Mini-Max principle.

- 1 Seat height adjustment / short perforated plate adapter
- 2 Seat height adjustment / long perforated plate adapter
- 3 Internal cross-brace
- 4 Neon Swing Away: with Swing-In and Swing-Out Function

Optimal adjustability

A wide variety of adjustment possibilities and options available through a minimum of handles enables the Neon to be adapted to the most varied requirements. The positioning of the back can be carried out while you are sitting in the wheelchair. In this way, you can feel every change immediately, so the user can find the most comfortable seating position for himself, without having to keep getting in and out, in the shortest possible time.



1



2



3

- 1+2 Undoing the adjusting screw
- 3 Adjusting the backrest

Back adjustment available as an option

QUICKIE

There is room for the Neon even in the smallest car.



Seat widths	from 30 – 50 cm
Seat depths	from 30 – 50 cm
Back heights	from 2.5 – 47.5 +/- 2.5 cm
27 individual frame colours	
6 individual sling colours	
2 frame angles	95°, 105°
2 frame offsets	0 cm, 2 cm
Footrests	Removable and swing-away footrests, on Neon Swing-Away
Weight	from 9.7 kg / 10.05 kg Neon Swing Away
Optional seat depth growth	2.5 and 5 cm (User weight max. 70 kg)
Suspension system	
Loading up to 42 cm seat width	115 kg/18 st
Loading from 42 – 50 cm seat widths	140kg/ 22 st

Now available – Lifestyle options:
26" lightweight wheels, 24" Offroad-wheels, 4" and 5" Design-castors;
24" and 26" Spinergy wheels

2005 Sunrise Medical
1005/2.000/000690277.UK/SD
Ref. A

STLPS
Sunrise Training and Education Programmes

KNOWLEDGE FOR THE FUTURE

**SUNRISE
MEDICAL**

Sunrise Medical Limited
High Street
Wollaston, Stourbridge
West Midlands DY8 4PS, England
Phone +44 (0) 1384 446622
Fax +44 (0) 1384 446644
www.sunrisemedical.com



**INNOVATE.
BUILD FOR THE FUTURE**