

XENON[®]

XENON[®] SA



We made the best even better. Constant improvements keep our product on the highest performance level.



// STRONGER ANGLE ADJUSTABLE BACK BRACKET

The material of the back bracket has been reinforced and the screws have been changed to stainless steel material.



// ADDITIONAL FRAME LENGTH

There is now a fourth frame length available for the Xenon SA, that reduces the total length of the wheelchair for certain seat depths. This also improves the lower leg position.



// IMPROVED FOLDABLE SEAT UPHOLSTERY

Now the tension of the seat upholstery is easier to adjust while mounted on the chair and the pocket sizes have been improved. Leading to a better folding



// HANGER STOP

The hanger stop for the Xenon SA prevents an unintentional opening of the knee lever wheel lock especially in combination with short seat depths.



// UPGRADED FOOTREST CLAMPING PIECE

for higher clamping forces.



// REDESIGNED SINGLE POST ARMREST RECEIVER

The connection between armrest and receiver is now free from backlash through two adjustment screws.