

# UltraLift



EASY TO LIFT, SIMPLE TO USE

Going out for the day in the car with your UltraLite scooter or powerchair couldn't be easier with the new UltraLift. The UltraLift from Electric Mobility is a simple mechanism that helps you lift your UltraLite into the boot without the worry of straining or having to ask for help.

When you want to go further with your UltraLite and remain independent, there is no better bootlift to choose than the UltraLift from Electric Mobility.


# UltraLift




## ULTRALITE™ VEHICLES

The UltraLift from Electric Mobility is easy to fit and simple to use. Easily attached to most tow bars\*, the UltraLift consists of a simple mechanism that uses power from the car's electrical system, cigarette lighter or the scooter's battery, eliminating the need for additional wiring. When you have stored your UltraLite, the lifting arm can quickly be taken apart and fits neatly in the boot. The UltraLift can be fitted for lifting from the left or right of a vehicle and has a lifting capacity of 110lb (50kg).


Take your independence further with the UltraLift from Electric Mobility.



Quick to dismantle and fits into your boot after use



Can be powered from the scooter battery or the car's cigarette lighter



Easy to use hand control



Full strapping kit is supplied



The UltraLift is easy to attach to most tow bars without the need for additional modifications\*

\* Tow bar is not included and must be installed by a competent fitter.

# EM BootLifts



When you need to lift a heavier scooter or powerchair, Electric Mobility gives you two different fixed bootlifts\*\* to choose from, both with a lifting capacity of 200lb (91kg).

### The EM 'stow away' BootLift

Ideal for occasional use, this fixed bootlift is wired directly into your car's electrics and can be easily dismantled into two parts for storage.

### The EM 'fold away' BootLift

Ideal for regular use, this fixed bootlift is both easily accessed and quick to fold away.



The EM 'stow away' BootLift dismantles completely – ideal for occasional use



The EM 'fold away' BootLift folds down as a complete unit – ideal for regular use

\*\* Please check with your local garage that your vehicle is suitable for an EM BootLift to be fitted.



For further details on the UltraLift and the EM BootLifts, or to discuss other models in the extensive Electric Mobility range –

Your local Electric Mobility dealer

## freephone 0800 252614

# ELECTRIC MOBILITY EURO

Electric Mobility Euro Limited, Canal Way, Ilminster, Somerset TA19 9DL  
Tel: 01460 258100 Fax: 01460 258125 email: sales@electricmobility.co.uk  
www.electricmobility.co.uk



Due to a policy of continual improvement, Electric Mobility Euro Limited reserve the right to change specifications without prior notice.