



**Bulletin No:** PU1230  
**Date:** 8 October 2018

## Heads up with Matrix Elan

**matrix**<sup>®</sup>  
SEATING SERIES

### New smaller head pad option

Following immense success of this highly adjustable headrest, Invacare has introduced an additional 4-point head pad. Now available in 2 sizes, for both adults and paediatrics.



Image: Featuring small vs. large 4-point head pad

Due to the aluminium plate construction, the Elan head pads can be formed around the individual's head. The popular 4-point head pad provides 4 points of support to the occipital and temporal areas of the head for optimum positioning, while carefully designed to fit around the ears so it doesn't feel intrusive or restrict hearing.

Small 4-Point size: 23Wx18Hcm

Large 4-Point size: 28Wx26Hcm

Please refer to your revised Matrix Elan prescription form for further details.


[www.invacare.co.uk](http://www.invacare.co.uk)



## Bulletin Communication

### Changes to points of adjustment

You will see the improved 2<sup>nd</sup> generation model transitioning through on your orders soon which features changes to the hardware joints. This change means that the unserviceable side has been capped off so that adjustment can only be made from the one side for smoother adjustments and ease of use.



The middle joint will now also feature security hex screws.




[www.invacare.co.uk](http://www.invacare.co.uk)




## Bulletin Communication

### Spare parts now available


Matrx Elan features 6 points of adjustment using an allen key. Two optional levers are provided with the Matrx Elan that allow for easy re-adjustment of the headrest or swingback function when hoisting or transferring. Following customer feedback of reissuing these high quality durable headrests in the future, Invacare have set up various spare parts to allow you to do this.



**Spare clamp kit**  
(code: 1629821)



**Spare knob kit**  
(code: 1629819)



**Spare screw kit**  
(code: 1629820)

For further details or a demonstration of Matrx Elan, please contact your local Seating Territory Manager .

[www.invacare.co.uk](http://www.invacare.co.uk)