

2017

TOPRO

**REHABILITATION
PRODUCTS**



TOPRO Limited

E-mail: info@topro.co.uk

Phone: 0844 739 1617

Fax: 0844 739 1615

www.topro.co.uk



2017

TOPRO REHABILITATION PRODUCTS give people with reduced mobility an easier everyday life

- Manufactured in Gjøvik, Norway by TOPRO
- Modern design to enhance and improve quality of life
- Developed together with healthcare professionals
- Manufactured with advanced technology
- Quality assurance according to NS-EN ISO 9001
- Environmental assurance according to NS-EN ISO 14001

TABLE OF CONTENTS

ODYSSEÉ 3
Travel rollator

TROJA 2G 7
Rollator

ROLLZ FLEX 11
Shopping rollator

OLYMPUS 13
Rollator

TROJA CLASSIC 17
Rollator

TROJA NEURO 21
Rollator

UMOTION 25
Rollator

TROJA WALKER 27
Forearm walker

TAURUS 31
Walker

ROLLZ MOTION 35
Rollator & wheelchair

**VERONA, MODENA
& SIENA** 37
Rise & recline chairs

BUTLER 43
Table

www.topro.co.uk

**TOPRO
ODYSSÉ**

TOPRO

New!



Made in Norway



THE TRAVEL ROLLATOR!

Safe, innovative two-way folding makes the TOPRO Odysseé perfect for transport and travel

CE

www.topro.co.uk



Instruction video:
<https://www.youtube.com/channel/UCeWbY5zdOQLZdoBM4TzovDw>



Stays securely upright
and takes little space when
folded one-way (sideways)



An indicator, in front
of the rollator, shows
when the unfolding
upwards is **safely in
the locked position**



Easy to fold with
a mechanism
for downwards
folding at a **user
friendly height**

Supportive travel companion

- Unique two-way folding (patent pending)
- Safe folding mechanism – cannot be folded downwards before completely folded sideways
- Compact when folded; 75 × 44.5 × 24.5 cm
- Practical transport bag is included
- Spacious, expandable shopping bag; 13–22.5 litres. Maximum weight capacity is 10 kilos. The fastening strap can be used to carry the shopping bag separately

Highly acclaimed features

- Different tyres are available which can easily be exchanged due to TOPRO Quick Release wheels
- Tilt function, edge guard, TOPRO Ergo Grips and a TOPRO Bell
- Safe, steady and lightweight
- Tested and approved according to EN ISO 11199-2 and biocompatibility tested according to EN ISO 10993-5
- Wide range of accessories are available

TOPRO ODYSSÉ is designed and approved for both indoor and outdoor use

7
Years
Warranty

TOPRO



Easy to transport in the specially designed bag
when folded two-way (sideways and downwards)



Small and handy
with the innovative
folding downwards

TOPRO ODYSSE

| | |
|---|---------------|
| Recommended user height | 150–200 cm |
| Maximum weight of user | 150 kg |
| Weight capacity shopping bag | 10 kg |
| Volume expandable shopping bag | 13–22,5 l |
| Handle height | 81.5–100 cm |
| Handles – Ergonomic (different holding positions – TOPRO Ergo Grip) | |
| Maximum length | 65 cm |
| Maximum width | 60 cm |
| Width between handles | 47 cm |
| Seat height | 62 cm |
| Diameter turning circle | 84 cm |
| Diameter wheels | 20 cm |
| Width wheels | 3.6 cm |
| Grip distance driving brake | 6.5 cm |
| Height when completely folded (sideways and downwards) | 44.5 cm |
| Height when only folded sideways | 81.5 cm |
| Width when folded | 24.5 cm |
| Length when completely folded (sideways and downwards) | 75 cm |
| Transport weight (excl. detachable wheels and shopping bag) | 6.0 kg |
| Total product weight (incl. detachable wheels and shopping bag) | 8.2 kg |
| TOPRO ART. NO. | 815000 |



Adjustable back support
with padding

814607



Tray

814728



Anti slip mat for tray

814059



Crutch holder

814457



Bottle holder

814043



Flower stickers
set of 4

814141



Nameplate

814024



Lock

814023



Lamp

814608



Umbrella
with attachment

814068



Off-road wheels for IBS,
set of 4

814649



Studded tyres for IBS,
set of 2 rear wheels

814670



Ice spurs for shoes, M (36-40)
L (shoe size 41-46)

814029

814035



One-handed brake

814026



Guiding handle
for e.g. visually impaired

814014



Rear bag with zipper

814046



THE NEW GENERATION OF ROLLATORS!





Increased safety, flexibility and the best possible comfort

The TOPRO Troja Classic revolutionised the rollator market in 2002. Since then the TOPRO product development team has developed a rollator based on research and feedback from users and therapists. The result is a 2nd generation TOPRO TROJA with several user friendly features and a strong focus on increased security, flexibility and comfort.

TOPRO TROJA 2G

7
Years
Warranty



INNOVATIVE!

Maximising user friendliness with various new functions

1 TOPRO IBS

A new developed interior braking system, which brakes on the rim instead of the surface of the wheel. This makes it possible to offer different tyres with a variety of tread patterns and softness, for different foundations.

2 TOPRO Quick Release

Easy to exchange or take off the wheels. Choose between:

3 Comfort Wheels Soft (PUR)

Comfortable and secure on unsurfaced roads, cobblestones etc.

4 Comfort Wheels Grip (TPE)

Manoeuvrable with good grip on smooth surfaces like flooring and tarmac (accessory).

5 TOPRO Off-road Wheels or Studded Tyres

For rugged terrain or icy roads (accessories).

6 TOPRO Ergo Grip

Unique ergonomic grip that makes it easier to stand from a sitting position. Provides a comfortable armrest and offers different holding positions when standing and resting.

7 TOPRO Bell

Attract attention when necessary. Easy to use without having to take your hands off the handles.

8 TOPRO Shopping Bag

Holds up to 10 kg and is detachable. It includes a shoulder strap for easy handling. The shopping bag has two pockets inside for small personal belongings and a lid to protect your items.

9 Reflective Panels

Be visible from all sides in the dark. The 3M Scotchlite Reflective panels increase safety.

In addition are the original features kept:

- Lightweight.
- Memory function for easy height adjustment of the handles.
- Tilt function and edge guard for safe and easy manoeuvring.
- Easy to fold.
- Stays firmly on its wheels when folded and takes little space.
- Tested and approved according to NS-EN ISO 11199-2.
- Biocompatibility tested according to EN ISO 10993-5.



Yasko Hupfauer
Physiotherapist at the University and rehabilitation clinic in Ulm, Germany

"Patients with reduced mobility and painful hands can hold the rollator with ease. The user walks WITHIN the rollator, not behind and the material of the grip is comfortable and provides good grip. By placing your hands on the top flat part of the Ergo Grip you will achieve an equilibrium reaction in the body."



Fabian Crüwell
Physiotherapist and therapy manager at Special Hospital Neresheim, Germany

"It gives a number of grip options. The grip requires more tension in the upper body, which is good for the patient. Rising up is easier, because you can support yourself better with this grip."



Nora Kern
NOI and Bobath instructor, Germany

"Sitting on the walker is more comfortable, since the arms can easily rest on the Ergo Grip."

TOPRO TROJA 2G

| | Medium | Small |
|--|-------------------|-----------------|
| Recommended user height | 150–200 cm | 135–170 cm |
| Maximum weight of user | 150 kg | 125 kg |
| Capacity shopping bag | 10 kg / 15 l | 10 kg / 15 l |
| Handle height | 78–100 cm | 67–86 cm |
| Handles – Ergonomic which offers different holding positions (TOPRO Ergo Grip) | | |
| Maximum length | 65 cm | 65 cm |
| Maximum width | 60 cm | 60 cm |
| Width between handles | 47 cm | 47 cm |
| Seat height | 62 cm | 54 cm |
| Diameter turning circle | 84 cm | 84 cm |
| Diameter wheels | 20 cm | 20 cm |
| Width wheels | 3.6 cm | 3.6 cm |
| Grip distance driving brake | 6.5 cm | 6.5 cm |
| Total dimensions L x W x H | 65 x 24 x 81.5 cm | 65 x 24 x 74 cm |
| Transport weight (excl. wheels and bag) | 6.0 kg | 5.6 kg |
| Total product weight (incl. wheels and bag) | 7.7 kg | 7.3 kg |



Turquoise Emerald

| | Turquoise Emerald | Pink Rose Sublime | Wine Red Tartan | White Baroque | Dark Grey | Silver Grey |
|--------|-------------------|-------------------|-----------------|---------------|-------------|-------------|
| Medium | 814600-1314 | 814600-1313 | 814600-1312 | 814600-1311 | 814600-1310 | |
| Small | 814601-1414 | 814601-1413 | 814601-1412 | 814601-1411 | | 814601-1400 |



TOPRO TROJA 2G comes with an elegant matching shopping bag and Comfort Wheels Soft Design on bags from the left; Rose Sublime, Grey, Baroque and Tartan



Adjustable back support with padding

814607



Tray

814728



Anti slip mat for tray

814059



Net basket detachable

814625



Shopping bag detachable

814622

Grey

814672

Rose Sublime

814642

Emerald

814643

Tartan

814641

Baroque



Rear bag with zipper

814046



Transport bag

814042



Crutch holder

814025



Bottle holder

814043



Flower stickers set of 4

814141



Nameplate

814024



Lock

814023



Umbrella with attachment

814068



Lamp

814608



Off-road wheels for IBS, set of 4

814649



Studded tyres for IBS, set of 2 rear wheels

814670



Ice spurs for shoes

M (36-40) 814029

L (41-46) 814035



One-handed brake

814026



Guiding handle for e.g. visually impaired

814014



Basket for oxygen bottle (only compressed, not liquid)

814009

New!



THE SHOPPER ROLLATOR!

Go shopping in comfort and style

Whether you are going to your local shop or planning a longer trip to the weekly market or supermarket the ergonomically designed Rollz Flex will ensure you can enjoy any shopping trip with confidence and in comfort. A unique feature of the Rollz Flex is the 25-litre shopping bag.



Denim Grey



Bright Purple

The large modern bag makes shopping easier than ever

The bag is available in two colours: Bright Purple and Denim Grey.

The Rollz Flex bag is made out of a sustainable and waterproof material and features a practical inner bag which enables you to easily unpack your shopping.



ROLLATOR ROLLZ FLEX

Rollz Flex features

- 1 Stable and fully adjustable**
The height adjustable handlebar has been designed to ensure you can walk steady and supported in the correct upright position by simply adjusting the height to meet your requirements. You can also turn the handlebar around to allow you to have the Rollz Flex either face forward or towards you depending on which is most comfortable.
- 2 Comfortable**
Comfortable ride and smooth suspension is provided by the cushioned tires.
- 3 Safe**
Innovative drum brakes for safe and balanced braking.
- 4 Compact**
The Rollz Flex is easily foldable for storage and transportation.
- 5 Manoeuvrability**
Very light steering and easy motion with the manoeuvrable front wheels.
- 6 Accessibility**
The threshold aid ensures the Rollz Flex can easily negotiate doorsteps and curbs.
- 7 Comfort**
The integrated seat and back support allows the user to rest at any time. Resting on a comfortable seating pillow and lean into the handlebar in forward position, which acts as back support.



ROLLZ FLEX

| | |
|-----------------------|---------------------------------|
| Volume shopping bag | 25 litres |
| Maximum load rollator | 125 kg |
| Minimum height | 160 cm |
| Maximum height | 190 cm |
| Dimensions unfolded | 60 × 62 × 80–100 cm (l × w × h) |
| Dimensions folded | 60 × 33 × 80–100 cm (l × w × h) |
| Seat width | 44 cm |
| Seat height | 62 cm |
| Diameter wheels | 20 cm |
| Weight rollator | 6.9 kg |

| | | |
|----------|---------------|------------|
| Colours | Bright Purple | Denim Grey |
| Art. No. | 3010RF0001 | 3010RF0002 |

Warranty:

5 years on the frame. 2 years on other parts with exception of normal wear and tear.



rollz
takes you further

**TOPRO
OLYMPOS**

TOPRO



THE OFFROADER!
A tough and practical rollator for an active life


Made in Norway

CE

www.topro.co.uk

**TOPRO
OLYMPUS**

7
Years
Warranty



Designed especially for outdoor use, but can also be used inside

- The tyres are made for comfortable driving on rough foundation, for example cobblestones or dirt roads.
- Detachable spacious bag with a cover
- Made from heavy duty materials which can tolerate all kinds of weather.
- Easy to clean.
- Tested and approved according to EN ISO 11199-2 (rollator standard) and EN ISO 10993-5 (biocompatibility).

Easy to unfold and a height adjustment with memory function

- Unique memory function sets the handles back to the appropriate height.

Safe to use with excellent brakes, kerb lifter and edge guard

- Safe to use and practical when crossing thresholds, kerbs etc.
- Anatomic handles.
- Easy to operate brakes.
- TOPRO IBS is a new developed interior braking system, which brakes on the rim instead of the surface of the wheel. This will make it possible to offer different tyres with a variety of tread patterns and softness, for different foundations.
- A variety of accessories that give extra safety in different situations.

Folds easily together and requires little space

- The rollator folds together by pulling the strap on the seat.
- Stands steady on its wheels when folded.
- Takes little space when being stored.
- Is practical and easy to transport.
- Excellent for users who travel by car, bus, train or airplane.

| Size | M medium | S small |
|-----------------------------|--------------|--------------|
| Colour | Dark Grey | Silver |
| Recommended user height | 150–195 cm | 140–170 cm |
| Maximum weight of user | 150 kg | 150 kg |
| Maximum capacity basket | 15 kg / 12 l | 15 kg / 10 l |
| Handle height | 76–95 cm | 71–87 cm |
| Handles | Anatomic | Anatomic |
| Maximum length | 76 cm | 76 cm |
| Maximum width | 70.5 cm | 70.5 cm |
| Width between handles | 50 cm | 50 cm |
| Seat height | 60 cm | 54 cm |
| Diameter turning circle | 94 cm | 94 cm |
| Diameter front wheels | 27.5 cm | 27.5 cm |
| Diameter rear wheels | 22.5 cm | 22.5 cm |
| Width on front/rear wheels | 4.2 cm | 4.2 cm |
| Grip distance driving brake | 6.5 cm | 6.5 cm |
| Height when folded up | 86 cm | 83 cm |
| Width when folded up | 27 cm | 29 cm |
| Length when folded up | 76 cm | 76 cm |
| Product weight | 8.9 kg | 8.8 kg |
| Art. No. | 814300 | 814305 |





814726



815212



814059



814046



814042



814040




814043



814608



814024



814023



814068




814661



814671



814026



814014



814009

**TOPRO
TROJA CLASSIC**

TOPRO



Made in Norway

THE ORIGINAL!

This rollator has been on the market for more than 12 years!



www.topro.co.uk

TOPRO TROJA CLASSIC

7
Years
Warranty



THE ORIGINAL DESIGN AND STILL VERY POPULAR!

Easy to fold, requiring little space when stored

- Pull the cord on the seat and the rollator folds together horizontally.
- Push the handlebars down into the sidebars.
- Unique memory function sets the handles back to the right height.
- Stands securely on its wheels when folded.
- Takes little space and is very practical to transport.
- Excellent for users who travel by car, train, bus or airplane.

Excellent brakes, kerb lifter and edge guard

- Easy to negotiate thresholds and kerbs etc.
- Safe to use.

Elegant and practical lightweight

- Plastic coated aluminium tube and plastic couplings.
- A total of 12 ball bearings in the wheels and forks provide excellent driving comfort.
- Can be hosed and cleaned.
- All wear and tear parts can be replaced.
- Brakes can be adjusted.
- Tested and approved according to EN ISO 11199-2 (rollator standard) and EN ISO 10993-5 (biocompatibility).
- TOPRO TROJA Classic is designed and approved for both indoor and outdoor use.



TOPRO TROJA CLASSIC

M

medium

S

small

X

extra small

| | | | |
|-----------------------------|----------------------|-----------------|-----------------|
| Recommended user height | 150–200 cm | 135–170 cm | 125–160 cm |
| Maximum weight of user | 165 kg | 125 kg | 80 kg |
| Maximum capacity basket | 10 kg / 20 l | 10 kg / 20 l | 5 kg / 11 l |
| Handle height | 78–100 cm | 67–86 cm | 63.5–79 cm |
| Handles | Anatomic | Anatomic | Anatomic |
| Maximum length | 65 cm | 65 cm | 65 cm |
| Maximum width | 60 cm | 60 cm | 52 cm |
| Width between handles | 47 cm | 47 cm | 37 cm |
| Seat height | 62 cm | 54 cm | 50 cm |
| Diameter turning circle | 84 cm | 84 cm | 79 cm |
| Diameter wheels | 20 cm | 20 cm | 20 cm |
| Width wheels | 3.6 cm | 3.6 cm | 3.6 cm |
| Grip distance driving brake | 6.5 cm | 6.5 cm | 6.5 cm |
| Total dimensions L × W × H | 65 × 24 × 80 cm | 65 × 24 × 74 cm | 65 × 24 × 65 cm |
| Weight | 6.8 kg | 6.6 kg | 6.1 kg |
| Area of use | Indoors and outdoors | | |

| Size \ Colour | Wine Red | White | Dark Grey | Silver |
|---------------|------------|------------|------------|------------|
| Medium | 814750-530 | 814750-540 | 814750-100 | 814750-120 |
| Small | | | | 814770-200 |
| Extra small | | | | 814740-000 |



TOPRO TROJA CLASSIC

Accessories can be ordered separately

ART. NO.



Adjustable back support with padding

814726



X
Tray

814728



X
Anti slip mat for tray

814059



X
Net basket detachable

814625



X
Shopping bag detachable

814622

Grey

814672

Rose Sublime

814642

Emerald

814643

Tartan

814641

Baroque




X
Rear bag with zipper

814046



Transport bag

814042



Crutch holder

814025



Bottle holder

814043



Flower stickers set of 4

814141



Nameplate

814024



Lock

814023



Umbrella with attachment

814068



Lamp

814608




Reflective panels

814624



Drag brake (pair)

814032



One-handed brake

814026



Guiding handle e.g. for visually impaired

814014



X
Basket for oxygen bottle (only compressed, not liquid)

814009



Ice spurs for shoes, size M(36-40) size L(41-46)

814029

814035



Forearm support (pair) for Troja Classic M Troja Classic S and X Only indoors on plain foundation

814829

814845



Stabilizer bar for forearm walker Troja Classic M and S Troja Classic X

814729

814730



Back support for forearm walker

814765

New!




Made in Norway

EXTRA SAFE WITH MORE CONTROL!

The TOPRO Troja – now available with a reverse braking system and laser light, designed and developed for people with Parkinson's and other neurological conditions



TOPRO TROJA NEURO

TOPRO TROJA NEURO Standard M



TOPRO TROJA NEURO Laser M



Extra safe

- The reverse braking system ensures the brakes are automatically activated
- The rollator will only move when the brake handles are pushed in
- If the brake handles are released the TOPRO Troja Neuro will stop immediately
- The reverse brake system will prevent falling
- Safe to sit down and stand up with the brakes automatically turned on
- Possible to deactivate the brakes completely by pushing the brake handles all the way down to the locked position. The rollator can then move freely
- Small turning circle makes it easy to use in small rooms like bathroom and kitchen
- The TOPRO Troja Neuro Laser model has in addition a green laser light that helps with "freezing" episodes which can occur with Parkinson's disease
- Easy to negotiate thresholds and kerbs etc. with the edge guard and tilt function
- Accessories for a variety of uses are available, e.g. one-handed brake and adjustable back support

Easy to fold, requiring little space when stored

- Pull the cord on the seat and the rollator folds together horizontally
- Push the handlebars down into the sidebars
- Unique memory function sets the handlebars back to the appropriate height
- Stays securely on its wheels when folded
- Takes little space and is very practical to transport
- Excellent for users who travel by car, train, bus or airplane

Practical and lightweight

- Plastic coated aluminium tube and plastic couplings
- Fabric consists of plastic and polyester
- A total of 12 ball bearings in the wheel and forks provide excellent driving comfort
- Can be steamed, hosed and washed, provided that the laser light unit and the switch box are removed

TOPRO TROJA NEURO is designed and approved for both indoor and outdoor use





TOPRO TROJA NEURO LASER



Laser light with a wireless switch

The TOPRO Troja Neuro Laser model comes with a laser light that can help a person who suffers from "freezing" episodes to move forwards. The laser light can easily be operated from the handle and the transfer is wireless to the laser light unit.



The laser light unit receives wireless signals from the switch on the handle and produces a green laser light which shows as a line between the rear wheels of the rollator. After 15 seconds the light is automatically turned off. If desired this can be adjusted. The green laser light is known to be six times more visible than red.

TOPRO TROJA NEURO

| | M medium | S small |
|--|--------------------|-------------------|
| Colour | Dark Grey | Silver |
| Recommended user height | 150–200 cm | 135–170 cm |
| Maximum weight of user | 165 kg | 125 kg |
| Maximum capacity basket | 10 kg / 16 l | 10 kg / 16 l |
| Handle height | 78–100 cm | 67–86 cm |
| Handles – Ergonomic which offers different holding positions (TOPRO Ergo Grip) | | |
| Maximum length | 65 cm | 65 cm |
| Maximum width | 60 cm | 60 cm |
| Distance between support points | 47 cm | 47 cm |
| Seat height | 62 cm | 54 cm |
| Diameter turning circle | 84 cm | 84 cm |
| Diameter front/rear wheels | 20 cm | 20 cm |
| Width front/rear wheels | 3.6 cm | 3.6 cm |
| Grip distance driving brake | 6.5 cm | 6.5 cm |
| Height when folded | 80 cm | 74 cm |
| Width when folded | 24 cm | 24 cm |
| Length when folded | 65 cm | 65 cm |
| Product weight STANDARD MODEL | 7.8 kg | 7.4 kg |
| Product weight LASER MODEL | 8.1 kg | 7.7 kg |

ART. NO. STANDARD MODEL

814746

814747

ART. NO. LASER MODEL

814743

814744



Adjustable back support with padding **814607**



Tray **814728**



Anti slip mat for tray **814059**



Rear bag with zipper **814046**



X Shopping bag detachable **814622** Grey **814672** Rose Sublime **814642** Emerald **814643** Tartan **814641** Baroque




Crutch holder **814025**




Lamp **814608**



Nameplate **814024**



Lock **814023**



Flower stickers set of 4 **814141**



Bottle holder **814043**




Umbrella with attachment **814068**



Transport bag **814042**



One-handed brake **814026**



Guiding handle for e.g. visually impaired **814014**



Basket for oxygen bottle (only compressed, not liquid) **814009**



*"Moving safely is
good medicine"
Prof. Dr. Bloem, neurologist*

UNIQUE STABILITY!

**When regular rollators do not meet the need, e.g.
with Parkinson's and other neurological conditions**

- Very safe and stable with reversed braking system.
- The rollator will only move when the brake handles are held in.
- Safe to sit down and get up as the brakes are always on.
- Single handed operating braking system.
- Small turning circle, convenient in small rooms like the bathroom.
- Integrated, adjustable drag brake for extra control.
- Green laser light to overcome "freezing" episodes.



VilleriusMedical



uMotion features

1 Reversed braking system

Prevents falling as the uMotion's brakes are always on and it will only move when the brake handles are held in. It stops immediately when the brake handles are released. The brakes can be operated with one hand; either brake handle works on both rear wheels.

2 Comfortable padded seat with back support

Sit safely and comfortable on the large seat with anti slip coating. Or use the seat to transport items.

3 Tension control

The desired resistance can be adjusted for extra control and support. The tension can also be adjusted asymmetrically.

4 Kerb support

The kerb support has three levels for high and low kerbs. The small wheels at the very end of the frame assist the uMotion over the kerbs and avoids tilting.

5 Stability

The U-shaped base provides unique stability and stimulates correct body posture.

6 Easy folding

Compact when folded for transportation or storage. Can be stored vertically.

7 Green laser light to overcome "freezing"

The freezing episode ceases when the user steps over the green laser line. Green laser light is six times more effective than red.



Accessories to be ordered separately
Shopping bag Art. No. 950025

UMOTION – rollator for special needs

| | |
|-------------------------|-----------------------------|
| Recommended user height | 150–185 cm |
| Maximum user weight | 125 kg |
| Handle height | 79–98 cm |
| Seat height | 54 cm |
| Seat width | 46 cm |
| Wheel diameter | 20 cm |
| Diameter turning circle | 85 cm |
| Width unfolded | 61 cm |
| Size folded | 95 × 63 × 46 cm (l × w × h) |
| Weight | 10 kg |

The uMotion is available in two variants; with or without laser light, see point 7. We do not recommend to mount the laser light on after purchase since it is a difficult task to undertake.

| | | |
|-----------------------------|----------|--------|
| uMotion with laser light | Art. No. | 950026 |
| uMotion without laser light | Art. No. | 950023 |

**TOPRO
TROJA WALKER**

TOPRO




Made in Norway

LIGHTWEIGHT WALKER!
Easy to fold, requiring little space when stored!

www.topro.co.uk

CE

TOPRO TROJA WALKER

7
Years
Warranty



ADJUSTABLE FOREARM SUPPORT!

Elegant and practical lightweight walkers

- Plastic coated aluminium tube and plastic couplings.
- A total of 12 ball bearings in the wheels and forks provide excellent driving comfort.
- Can be hosed and cleaned.
- All wear and tear parts can be replaced and brakes can be adjusted.
- Tested and approved according to EN ISO 11199-3 (walker standard).
- TOPRO TROJA FOREARM WALKER is designed and approved to be used indoors on flat surface.

Easy to fold, requiring little space when stored

- Pull the cord on the seat and the walker folds together horizontally.
- Push the handle bars down into the side bars.
- Unique memory function sets the handles back to the desired height.
- Stands securely on its wheels when folded.
- Takes little space and is very practical to transport.

Excellent brakes, kerb lifter and edge guard

- Easy to negotiate thresholds and kerbs etc.
- Safe to use.



TOPRO TROJA WALKER

M

medium

S

small

X

extra small

| | | | |
|-------------------------------|--------------------------|-----------------|-----------------|
| Recommended user height | 150–185 cm | 130–170 cm | 130–160 cm |
| Maximum weight of user | 135 kg | 125 kg | 80 kg |
| Maximum capacity shopping bag | 5 kg / 20 l | 5 kg / 20 l | 5 kg / 11 l |
| Height of support points | 98–117 cm | 84–100 cm | 84–96 cm |
| Diameter handles | 3.5 cm | 3.5 cm | 3.5 cm |
| Maximum length | 65 cm | 65 cm | 65 cm |
| Maximum width | 63 cm | 63 cm | 53 cm |
| Width between support points | 45 cm | 45 cm | 35 cm |
| Seat height | 62 cm | 54 cm | 50 cm |
| Diameter turning circle | 84 cm | 84 cm | 79 cm |
| Diameter wheels | 20 cm | 20 cm | 20 cm |
| Width wheels | 3.6 cm | 3.6 cm | 3.6 cm |
| Grip distance driving brake | 6.5 cm | 6.5 cm | 6.5 cm |
| Total dimensions L × W × H | 65 × 33 × 102 cm | 65 × 33 × 88 cm | 65 × 33 × 88 cm |
| Weight | 10.1 kg | 9.8 kg | 9.3 kg |
| Area of use | Indoors on flat surfaces | | |
| Art. No. | 814749 | 814742 | 814741 |



TOPRO TROJA WALKER

Accessories can be ordered separately

ART. NO.



Back support **814765**



X Tray **814728**



X Anti slip mat for tray **814059**



X Net basket detachable **814625**



X Shopping bag detachable **814622** Grey **814672** Rose Sublime **814642** Emerald **814643** Tartan **814641** Baroque



X Rear bag with zipper **814046**



Crutch holder **814025**



Bottle holder **814043**



Flower stickers set of 4 **814141**



Nameplate **814024**




Lock **814023**



Drag brake (pair) **814032**



Stabilizer bar for forearm walker Troja Classic M and S **814729**
Troja Classic X **814730**



Guiding handle **814014**
e.g. for visually impaired



X Basket for oxygen bottle (only compressed, not liquid) **814009**

**TOPRO
TAURUS**

TOPRO



SAFE AND STEADY!

Ideal walkers for domestic and care environments


Made in Norway

CE

www.topro.co.uk

TOPRO TAURUS

Steady and safe walkers for domestic and care environments

- Easy to read height scale and easy to adjust the height.
- High enough for tall users (up to 2.10 m).
- Tested and approved according to ISO 11199-3 (walker standard) for use indoors on flat surfaces.
- The TOPRO TAURUS E-models are tested and approved according to IEC 60601-1 (electrical safety) and IEC 60601-1-2 (electromagnetic compatibility) and EN ISO 10993-5 (biocompatibility).

Maintenance free material

- Aluminium with polyester surface and polyurethane.
- Easy to clean.
- Non-corrosive.
- All wear and tear parts can be exchanged.

Two height adjustment methods:

TOPRO TAURUS H – Hydraulic height adjustment.

TOPRO TAURUS E – Electrical height adjustment.

Two levels of equipment:

BASIC models have standard equipment with adjustable arm pads and foot brakes on all wheels. Delivered partly assembled.

PREMIUM models have in addition



| Model | Art. No. |
|------------------------|----------|
| TOPRO TAURUS H Basic | 814785 |
| TOPRO TAURUS H Premium | 814787 |
| TOPRO TAURUS E Basic | 814790 |
| TOPRO TAURUS E Premium | 814789 |

TOPRO TAURUS H

- Gas springs simplify the height adjustment.
- Can be used in the shower.

Product weight 18 kg

Maximum weight of user 150 kg

Recommended user height 140–210 cm



Main image: TOPRO TAURUS H Basic.



TOPRO TAURUS E

- Electric motor makes the height adjustment very user friendly.
- Can also be used as raising support.

Product weight 22 kg

Maximum weight of user 150 kg

Recommended user height 140–210 cm



Foot plate

814645



Forearm side support (pair)

814647



Brake handle left side (instead of right) only for Premium models

814846



Adjustable handles

814001



Tray

814810



Height adjustment panel (Concens)

814646



Wall mount for the Taurus battery charger (Concens)

814085



Small basket

814821



Infusion stand

814793



Stand for body fluid-boxes

814034



Oxygen bottle holder

814820



EASY TO SWITCH!

Rollator and wheelchair in one provides you more freedom

The Rollz Motion is an ergonomically designed rollator with a strong but lightweight frame, which enables you to walk steadily, upright and comfortably. If you want to rest along the way, the Rollz Motion converts from a rollator to a wheelchair at a flick of the wrist, and of course, just as easily back again. The Rollz Motion is developed and produced by the Dutch company Rollz.



reddot design award
honourable mention 2012



product
design award

2012

The Rollz Motion is available in four colours:
Purple (limited edition), Orange, Ice Blue and Light Grey.



ROLLATOR ROLLZ MOTION

From rollator to wheelchair in an instant!

- Before you leave home, click the Rollz wheelchair attachment into place on the rollator.
- Pull the backrest into position in one flowing motion.
- Turn the handles from the rollator position into the push position, and adjust them to the appropriate height for the person pushing the wheelchair.
- Click the footrests into place at the front.

Continue your journey in comfort and style!



Rollz Motion features

1 Stability

Walking upright, steadily and supported between the ergonomic hand grips that are adjustable in height.

2 Safe

Innovative drum brakes for safe and balanced braking.

3 Ease

Cross doorsteps and curbs easily with the threshold aid.

4 Manoeuvrable

Light steering and easy motion with the manoeuvrable front wheels.

5 Comfortable

Comfortable ride and natural suspension thanks to the foam-filled rear tyres.

6 Convenient

Resting on a comfortable cushioned seat.

7 Compact

The Rollz Motion is easily foldable into a compact package.



Accessories can be ordered separately

| | Art. No. |
|-------------------------------|------------|
| A Shopping bag | 1020RM0004 |
| B Cane holder | 1020RM0003 |
| C Cup holder | 1020RM0002 |
| D Travel bag | 1020RM0010 |
| E Footrests size S | 1030RM0008 |
| Footrests size L | 1030RM0010 |
| Footrests size M are included | |
| Reflection stickers | 1020RM0005 |
| Backrest | 1020RM0012 |

ROLLZ MOTION – rollator and wheelchair in one!

| | |
|-----------------------|--|
| Volume shopping bag | 25 litres |
| Maximum load rollator | 125 kg |
| Minimum height | 160 cm |
| Maximum height | 190 cm |
| Dimensions unfolded | 60 × 62 × 80–100 cm (l × w × h) |
| Dimensions folded | 60 × 33 × 80–100 cm (l × w × h) |
| Seat width | 44 cm |
| Seat height | 62 cm |
| Diameter wheels | 20 cm |
| Weight rollator | 6.9 kg |
| Colours | Purple Orange Ice Blue Light Grey |
| Art. No. | 2010RM0006 2010RM0001 2010RM0002 2010RM0004 |

Warranty: 5 years on the frame.

2 years on other parts with exception of normal wear and tear.



RELAX IN COMFORT!

Sit, rest or watch TV comfortably in a TOPRO rise and recline chair

It provides comfortable and safe support for your back and legs with the ergonomic seating and backrest.

With a built-in stand-up assistance, you can adjust the level of support you need to stand up or sit down with confidence.

- The remote control buttons are large and easy to use.
- The chairs are easy to clean and maintain because of the microfibre material or faux leather.
- All models in the microfibre material are available in a cream, caramel or mocca colour and the faux leather is available in brown, cream or red.
- The chair comes in two parts making it easy to assemble.

All Chairs and accessories are available in microfibre and faux leather in a variety of colours



VERONA

Our most versatile chair. With a dual motor for a more flexible adjustment and very practical

MODENA

Elegant, classy, stylish and modern

SIENA

Safe and comfortable with a wing shaped back offering good lateral support

| | | | | |
|--|--------------------------------|--------------------|--------------------|--------------------|
| | Recommended user height | Approx. 165–185 cm | Approx. 160–180 cm | Approx. 160–180 cm |
| | Recommended max. user weight | 150 kg | 150 kg | 150 kg |
| | Chair height | 107 cm | 105 cm | 103 cm |
| | Chair width | 77 cm | 83 cm | 80 cm |
| | Chair length, resting position | 175 cm | 162 cm | 162 cm |
| | Back height | 59 cm | 61.5 cm | 57.5 cm |
| | Seat height | 48 cm | 43.5 cm | 45.5 cm |
| | Seat width | 52 cm | 51 cm | 51 cm |
| | Seat depth | 52 cm | 52 cm | 54 cm |
| | Back tilt | Approx. 180° | Approx. 158° | Approx. 158° |

Intuitive remote control VERONA

- Sit
- Relax
- Get up
- Footrest – can be adjusted separately
- Backrest – can be adjusted separately
- Verona is adjustable to a horizontal sleeping position



Intuitive remote control MODENA and SIENA

- Two options for two different positions

- Child protection pin



VERONA

The most comfortable and adjustable chair from TOPRO


- Our top of the range model with dual motor provides control of the backrest and the footrest independently.
- The dual motor gives you reliable and safe support when sitting down and standing up.
- The remote control can be mounted on the left or the right hand side of the chair.
- A comfortable neck pillow made of memory foam is included as standard.
- With four wheels, the chair is easy to move without having to lift.
- Enables a relaxing posture when sitting with outstretched legs and the back straight.
- The seat is adjustable to a horizontal sleeping position.
- Parking brake on the rear wheels.
- Accessories (available separately)
 - A memory foam lumbar cushion for the lower back.
 - Arm rest covers (pair).



| The VERONA is available in microfibre and faux leather in a variety of colours | Microfibre | | Art. No. | Faux leather | | Art. No. |
|--|------------|---------|----------|--------------|---------|----------|
| | Mocca | Caramel | Sand | Brown | Cream | Red |
| VERONA chair Neck support pillow included | 8734-MB | 8734-DB | 8734-SA | 8734-BR | 8734-CR | 8734-RE |
| Lumbar cushion | 0039-MB | 0039-DB | 0039-SA | 0039-BR | 0039-CR | 0039-RE |
| Armrest covers | 0037-MB | 0037-DB | 0037-SA | 0037-BR | 0037-CR | 0037-RE |

Accessories

The lumbar cushion and armrest covers are available in all colours and materials



For measurements and other product details, please see page 38.

MODENA

Modern, safe and comfortable

- The Modena provides perfect support to help one stand up or sit down.
- The Modena is fitted with a motor to help you safely transfer from the sitting to standing position.
- The remote control can be mounted on the left or right hand side and is equipped with a child safety lock pin.
- The two wheels make the chair easy to move.
- Accessories (available separately)
 - A memory foam lumbar cushion for the lower back.
 - Supportive neck pillow made of memory foam.
 - Arm rest covers (pair).



| The MODENA is available in microfibre and faux leather in a variety of colours | Microfibre | | Art. No. | Faux leather | | Art. No. |
|--|------------|---------|----------|--------------|---------|----------|
| | Mocca | Caramel | Sand | Brown | Cream | Red |
| MODENA chair | 8432-MB | 8432-DB | 8432-SA | 8432-BR | 8432-CR | 8432-RE |
| Neck support pillow | 0040-MB | 0040-DB | 0040-SA | 0040-BR | 0040-CR | 0040-RE |
| Lumbar cushion | 0039-MB | 0039-DB | 0039-SA | 0039-BR | 0039-CR | 0039-RE |
| Armrest covers | 0038-MB | 0038-DB | 0038-SA | 0038-BR | 0038-CR | 0038-RE |

Accessories

The neck support pillow, lumbar cushion and armrest covers are available in all colours and materials



For measurements and other product details, please see page 38.

SIENA

Safe and comfortable

- The wings give the chair a softer look and provide more security when leaning.
- The motor controls the footrest and backrest , making it easy to manoeuvre.
- The remote control can be mounted on the left or on the right hand side and is equipped with a child safety lock pin.
- The Siena is fitted with a motor to help you safely transfer from the sitting to standing position.
- The two wheels make the chair easy to move.
- Accessories (available separately)
 - A memory foam lumbar cushion for the lower back.
 - Supportive neck pillow made of memory foam.




| The SIENA is available in microfibre and faux leather in a variety of colours | Microfibre | | Art. No. | Faux leather | | Art. No. |
|---|------------|---------|----------|--------------|---------|----------|
| | Mocca | Caramel | Sand | Brown | Cream | Red |
| SIENA chair | 8431-MB | 8431-DB | 8431-SA | 8431-BR | 8431-CR | 8431-RE |

| | | | | | | |
|---------------------|---------|---------|---------|---------|---------|---------|
| Neck support pillow | 0040-MB | 0040-DB | 0040-SA | 0040-BR | 0040-CR | 0040-RE |
| Lumbar cushion | 0039-MB | 0039-DB | 0039-SA | 0039-BR | 0039-CR | 0039-RE |

Accessories

The neck support pillow and lumbar cushion are available in all colours and materials



For measurements and other product details, please see page 38.





A great help around the house

The TOPRO Butler provides support when walking, as well as being a serving station or extra table. It is practical and safe when moving items around the home, e.g. carrying food and drinks between the kitchen and the living room or carrying laundry from the washroom.

- Quick to assemble.
- Easy height adjustment.
- Turning wheels with brakes.
- Made of polyester coated aluminium.
- Table top is made from hard durable plastic.
- Easy to clean.

| | |
|--------------------------|------------|
| Width | 54 cm |
| Depth | 42 cm |
| Handle height | 71–102 cm |
| Height of table top | 49–80 cm |
| Table top, width × depth | 47 × 36 cm |
| Diameter wheels | 8 cm |
| Width wheels | 2.5 cm |
| Weight | 4.2 kg |
| Maximum load | 7 kg |

| | |
|--------------|----------|
| | Art. No. |
| TOPRO BUTLER | 814041 |

| | |
|---------------------------------------|--------|
| Accessories, to be ordered separately | |
| Basket | 814530 |



2017

TOPRO

**REHABILITATION
PRODUCTS**



TOPRO Limited

E-mail: info@topro.co.uk

Phone: 0844 739 1617

Fax: 0844 739 1615

www.topro.co.uk