

Invacare® | ***Flo-tech® Contour Lo-Back
Cushion***



Optimised surface support

The *Invacare Flo-tech Contour Lo-back* cushion has been specifically designed to have a lower rear edge allowing the client to sit slightly reclined. This gives postural support whilst reducing sacral area interface pressures and without additional and unwanted forces pushing forwards. This cushion is suitable for those 'At Risk'* of developing pressure ulcers.



Yes, you can.®



Flo-tech Contour
Lo-Back



Features and Options

Enhanced postural position

The **Lo-back cushion** reduces the need for a wedged cushion, lumbar roll or reclined backrest and minimises the tendency to slide forward. This reduces shear and friction forces and promotes an enhanced postural position.

The cushion is most effective when used in conjunction with an angle-adjustable back support and can also be used with a reclining back or tilt in space frame for those with extra postural needs.



Ingress Resistant Cover

Cover is multi stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Range of cushions covers available – Black Smooth, Towelling and PU covers.



High quality resilient foam offering high levels of pressure reduction and product longevity.

Technical data



Flo-tech Contour Lo-Back Cushion A range of 83 sizes



Height: 103 mm



Weight: 1.38 kg



Max. User Weight: 152 kg (24 stone)



Cleaning: 80°

Water Resistant
Vapour Permeable

Manufactured to European Standards:
BS EN 597-1 Cig, BS EN 597-2 Match, and BS 7176 (Crib 5)

Colours



Black spacer (Option)



Black Smooth (Optional)



Black Wipedown (Option)

INVACARE® Limited
Pencoed Technology Park - Pencoed - Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Flo-tech Contour Lo-Back Cushion
UK - 01/2011 - 1505626

* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

