

Invacare® | **Flo-tech® Image**



3 Year warranty

Slimline contoured cushion with gel sac

The *Invacare Flo-tech Image* cushion uses a combination of a slim, contoured foam with a large gel sac covering the sacral area. This increases even weight distribution and immersion of the Ischial Tuberosities. The design delivers a high degree of protection to clients who are considered to be at 'High Risk'* of developing pressure ulcers.



Yes, you can.®

Slimline contoured foam and gel combination cushion

Some clients prefer a slimline cushion, finding it more compatible with their daily living activities. The *Invacare Flo-tech Image* combines the benefit of a contoured foam upper surface with a large gel sac covering the ischial and sacral area, offering increased support and reducing the occurrence of peak pressures.



Contoured seating surface with gel sac

Offers improved stability and postural support. The moulded foam castellations and gel combine effectively to reduce shear and friction forces on tissue, and help to prevent the client sliding forward in the seat.

Features and Options



Water Resistant Cover

Cover is multi stretch, water-resistant and vapour-permeable. Cover features an anti-slip base and carry handle.

Available in a wide range of sizes

To view the cushion size range, please visit www.invacare.co.uk.

Technical data



A range of 50 sizes



Height: 7.8 cm / 3"



Weight: 1.82 kg



Max. User Weight: 140 kg (22 stone)



Cleaning: 80°

Flo-tech Image

Water Resistant
Vapour Permeable

Manufactured to European Standards:
BS EN 597-1 Cig, BS EN 597-2 Match, and BS 7176 (Crib 5)

Colours



Black Smooth



Black Wipedown



Black Towelling



Black Spacer Fabric



Mulberry

Warranty 3 Years

INVACARE® Limited
Pencoed Technology Park - Pencoed - Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Flo-tech Image - UK - 10/2010
1505623



Yes, you can.®

* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

