



INCREASED
CLINICAL FUNCTIONALITY

Invacare[®] Rea[®] Flex 3



Confidence in posture solutions

Trusted Solutions, Passionate People



Yes, you can.[®]

FLEX 3

Discover the Flex 3 system

The Flex 3 system has been specially designed to **fit a variety of people**, allowing a seat to be **tailored** to different body shapes and disabilities. Flex 3 aims to **improve posture and stability** which in turn increases activity levels and quality of life.



Compatibility between wheelchairs

With the same attachment on all Invacare passive wheelchairs*, the Flex 3 system provides the **flexibility** to combine backrest cushions with backrest types across its range. This allows for **backward compatibility** and **enhances refurbishment**.



* Excluding Clematis® and Azalea® Minor

Find the right wheelchair

Rea Azalea® Range

The **Rea Azalea** range is the natural choice for the moderately active client as well as for the more passive and offers a **25° tilt** with a **30° low shear recline** function. The range offers a solution for a wide variety of requirements.

Azalea	Adapts to a wide range of seating and positioning needs.
Azalea Base	Accommodates alternative seating solutions.
Azalea Minor	Smaller dimensions for teenagers and petite adults.
Azalea Tall	Higher seat height and longer seat depths for taller people.
Azalea Max	Meets the specific needs of the larger person.



Rea Dahlia® Range

The **Rea Dahlia** is manoeuvrable, slim and adjustable with impressive tilt options and a **30° low shear recline** function. This versatile wheelchair is perfect for those seeking mobility with comfort and pressure relief.

Dahlia 45 - Impressive **45° tilt** delivers excellent pressure relieving properties and supports the most vulnerable areas. To obtain pressure relief using only the tilt function, a wheelchair must go to the full 45° angle.

Dahlia 30 - **30° tilt** version with lower floor to seat height for people who can readily propel themselves. With a self-tilt function available, clients can actively change posture themselves when required.



Now select the backrest system

Easy backrest guide



Cover Lateral
The Cover Lateral follows the contoured shape set by the tension adjustable straps to give the client a good level of support.



Laguna
The Laguna backrest cushion offers moderate lateral support and gives the client stability and support.



Mistral 2
The Mistral 2 backrest cushion with its deep, reinforced wings and a lower, shaped upper section offers the client a distinct waist support and assists shoulder movement.



Passad 2
With its high shoulder support, the Passad 2 offers support for client with a weak posture. The high wings are reinforced.



Shoulder High
The Shoulder High backrest cushion is the same as Passad 2 but is 100 mm higher, and has been designed specially for taller clients.

Proposed backrest options

Plate



Short

Plate



Tall

Either plate or tension adj



Average

Plate



Wide

Tension adj



Wide above

Tension adj



Wide below

Tension adj



Larger in front

Tension adj



Larger behind

Tension adj



Kyphosis back

Tension adj



S-shaped back



Backrest Tension Adjustable

The backrest can be tailored to an individual shape by adjusting the straps. It is angle, width and height adjustable to provide a high level of support and comfort.



Backrest Plate

A width and height adjustable, firm and stable support plate.



Take a look at our Flex 3
video on our Invacare Europe
YouTube channel called
Invacare Flex 3 Seating System.



**If further postural or seating support is needed, we
have a wide range of accessories available. Take a look
at our Rea accessory brochure for more information.**



Invacare Limited
Pencoed Technology Park - Pencoed - Bridgend
CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed
to be correct at time of print.
Invacare reserves the right to alter product
specifications without prior consultation.
Flex 3 - UK - 11/2016 - 1553824

Trusted Solutions,
Passionate People



Yes, you can.