

Product description - All included
<p>Each order is composed of</p> <ul style="list-style-type: none"> • 1 power add-on • 1 small joystick • 1 anti-tip • 1 Li-ion battery • 1 charger • 1 installation kit <p>24 months warranty</p>

Dealer Information
Date of order: / /
Purchase Order No.
Dealer:

POWER ADD-ON	Benoit Light Drive 3 Å5,295.00 rrp			
COLOUR	● Black			
LENGTH OF TUBE ⁽²⁾	345 mm <i>between 42.5 to 48.5 cm track width ⁽¹⁾</i>	410 mm <i>between 49 to 55 cm track width ⁽¹⁾</i>	475 mm <i>between 55.5 to 61.5 cm track width ⁽¹⁾</i>	540 mm <i>between 62 to 68 cm track width ⁽¹⁾</i>
Joystick side	Left		Right	
WHEELCHAIR	Seat width:		Wheelbase width: <small>Centre to centre of tyre</small>	
	Make:		Model:	
	Type of backrest:	Fixed Folding Reclining	Type of wheelchair:	Folding frame Rigid frame
CHOICE OF BATTERY	484 Wh Li-ion battery		282 Wh Li-ion battery → <i>Not recommended with XL and Light Drive Plus</i>	

(1) Distance between the centers of the left and right rear wheels. (2) Please refer to the table «Wheelchair tube compatibility». → Comply with applicable IATA regulations

OPTIONS					
OPTION	ITEM NUMBER	DESIGNATION	PHOTO	QUANTITY	PRICE
JOYSTICK AND CONTROL MODULE	MAC-Z03	Additional controller bracket + fixation set			£97.78
	MAC-Z09	Dual attendant control			£563.00
	MAC-Z08	Additional controller bracket for dual attendant control			£129.00
ACCESSORIES	500-068	Car charger for Lithium-ion battery (cigarette lighter socket)			£261.05
	400-095	Mudguards			£95.69 (pair)
	750-010	Schwalbe Marathon Plus tyres 24"x1			£150.50 (pair)
	MAB-Z02	282 Wh Lithium-ion battery →			£1,044.14
	MAB-Z07	484 Wh Lithium-ion battery			£1,392.19
Susurro Installation	LD3-INSTALL	Benoit Installation carried out by our technician at Susurro facility. - Dealer responsible for shipping chair to OX29 8LJ - Susurro responsible for palletised return shipping of completed project			£295.00
Notes:					

Please send orders to orders@susurro.co.uk

WHEELCHAIR TUBE COMPATIBILITY


Update : Jan-2026

Table for tube selection according to **TRACK WIDTH⁽¹⁾**

Tube length (mm)	tube 345	tube 410	tube 475	tube 540	Custom-length tube	
TRACK WIDTH (cm)	42,5 43 43,5 44 44,5 45 45,5 46 46,5 47 47,5 48 48,5 49 49,5 50 50,5 51 51,5 52 52,5 53 53,5 54 54,5 55 55,5 56 56,5 57 57,5 58 58,5 59 59,5 60 60,5 61 61,5 62 62,5 63 63,5 64 64,5 65 65,5 66 66,5 67 67,5 68 68,5 69 69,5 70 70,5 71 71,5 72 72,5 73 73,5 74 74,5 75 75,5 76 76,5 77 77,5 78 78,5 79 79,5 80 80,5 81 81,5 82 82,5 83 83,5 84 84,5 85 85,5 86 86,5 87 87,5	Size not suitable for XL				

Table for tube selection according to **SEAT WIDTH**

Tube length (mm)	inertitude ⁽²⁾	tube 345	tube 345 or 410	tube 410	tube 410 or 475	tube 475	tube 475 or 540	tube 540	Tube 540 or custom-length tube	Custom-length tube	
SEAT WIDTH (cm)	30 30,5 31 31,5 32 32,5 33 33,5 34 34,5 35 35,5 36 36,5 37 37,5 38 38,5 39 39,5 40 40,5 41 41,5 42 42,5 43 43,5 44 44,5 45 45,5 46 46,5 47 47,5 48 48,5 49 49,5 50 50,5 51 51,5 52 52,5 53 53,5 54 54,5 55 55,5 56 56,5 57 57,5 58 58,5 59 59,5 60 60,5 61 61,5 62 62,5 63 63,5 64 64,5 65 65,5 66 66,5 67 67,5 68 68,5 69 69,5 70 70,5 71 71,5 72 72,5 73 73,5 74 74,5 75 75,5 76 76,5 77	Size not suitable for XL		May not be compatible with XL >Check with BENOIT SYSTEMES							



(1) Track width : Distance between the centers of the left and right rear wheels (approximately at the installation height of the power add-on)

(2) May not be compatible : you can spread the wheels and recheck the track width