

# Invacare Matrx Flo-tech Solution

The *Invacare* **Matrx Flo-tech Solution** cushion uses a combination of a slim, shaped foam cushion base and a dual layer of gel sacs over the ischial and sacral area. The dual gel sacs and the unique moulded foam shape of the solution increases the support surface area. This delivers a superior level of protection for those patients identified as being at Very High Risk of developing pressure ulcers, or for those with already established tissue damage.

# Slimline contoured foam base with large dual gel sac covering the entire top surface of the cushion

The contact area between the user and the foam cushion is protected, as the gel re-distributes pressure away from bony prominences. The additional layer of gel is incorporated under the ischial and sacral area to further enhance the pressure reducing qualities of this cushion, making this cushion suitable for Very High Risk users.

# Slimline contoured cushion

Some users prefer a slimline cushion, finding it more compatible with their daily living activities. The low profile is

particularly beneficial for the facilitation of side transfers.

### **Features and Options**



**Ingress Resistant Cover** Cover is two-way stretch, waterresistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Available in 44 sizes The *Invacare Matrx Flo-tech Solution* is available in in a range of 44 sizes, from 38 x 41 to 51 x 46 cm (380 x 410 to 510 x 460 mm).



# Invacare Matrx Flo-tech Deep Solution option

Option of a 10 cm (100 mm) "Deep" foam base is also available in 12 sizes.



#### Ultra-Fresh<sup>™</sup> added to all *Flotech* foam for antimicrobial and odour protection

Antimicrobial treatment that has been added to all *Flo-tech* foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.