

XT4



03/2025


ORDER FORM

Errors and omissions excepted..

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town/PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Art. No.:		
Purchase Order No.:			User weight: kg		
Reimbursement Code:			Base Price		
Order date: Day Month Year			£ 16,950		

Client Measurements

a) Top Of Head To Seat	mm
b) Top Of Shoulder To Seat	mm
c) Underarm To Seat	mm
d) Elbow To Seat	mm
e) Top Of Pelvis To Seat	mm
f) Shoulder Width	mm
g) Trunk Width	mm
i) Hip Width	mm
j) Actual Sitting Depth	mm
m) Foot Length	mm
n) Lower Leg Length	mm
o) Actual Seat Width	mm



Technical Data

User Weight Capacity	160 kg with seat lift and tilt	Ground Clearance	125 mm
Base Width	700 mm	Drive Wheels & Castors	Off-road pneumatic 145 x 70, black rim
Chassis length	1000 mm	Charger	220 V/12 A
Drive Control	R-Net (120A Power Module)	Battery Size	115 Ah
Drive Motors	10 km/h, 4 x 700 W, 4-Pole Motors	Turning radius	1410 mm
Standard configuration	Including LED lights, Crash Tested Vehicle Tie-Downs, Maintenance-free 115 Ah batteries, 12 A charger, MPS back system, MPS seat system		
Tyres	Maximum tyre pressure depending on user weight from 0.25 to 0.55 bar (see table in instruction manual)		

Drive Unit

Power Wheelchair

MAG010013	XT4	£
		16,950

Colour Selection

MAG6094	XT4 Midnight Black	£
MAG6095	XT4 Sunset Red	0
MAG6096	XT4 Horizon Blue	0

Puncture Resistance Options

MAG7050	Pneumatic	£
MAG7056	Extreme XT4 Puncture Resistant Kit	0
MAG7058	Neverflat for XT4	600
		160

Batteries

MAG6285	Battery, 115 Ah, Maintenance free	£
		0

Chargers

MAG160105	Charger, Automatic, 220 V / 12 A, UK	£
		0

Seat Tilt Modules

MAG3502-KIT	Fixed Frame (only lighting)	£
MAG5335	Fixed Frame (with Power Options)	0
MAG5127	Tilt (0 - 50°)	0
MAG5126	Lift and Tilt (300 mm lift, 0 - 45°)	2950

Seat Height

MAG6552	Seat Height: Leave as is Max Lower Leg Length 370 mm	£
MAG6553	Seat Height: +25mm front, +25mm rear Added Seat Height Versus 480 mm Standard Seat Height	0
MAG6554	Seat Height: +50mm front, +50mm rear Added Seat Height Versus 480 mm Standard Seat Height	0

Seating Systems

MAGREHABSEATPAN	Rehab Seat Pan	£
		0

Rehab Seat Sizes

Seat width Depth minus Width Must Not Exceed 50 mm

MAGSW305MM	305 mm	MAGSW430MM	430 mm	MAGSW560MM	560 mm
MAGSW330MM	330 mm	MAGSW455MM	455 mm	MAGSW585MM	585 mm
MAGSW355MM	355 mm	MAGSW480MM	480 mm	MAGSW610MM	610 mm
MAGSW380MM	380 mm	MAGSW510MM	510 mm		
MAGSW405MM	405 mm	MAGSW535MM	535 mm		

Seat depth Depth minus Width Must Not Exceed 50 mm

MAGSD305MM	305 mm	MAGSD430MM	430 mm	MAGSD560MM	560 mm
MAGSD330MM	330 mm	MAGSD455MM	455 mm	MAGSD585MM	585 mm
MAGSD355MM	355 mm	MAGSD480MM	480 mm	MAGSD610MM	610 mm
MAGSD380MM	380 mm	MAGSD510MM	510 mm		
MAGSD405MM	405 mm	MAGSD535MM	535 mm		

Seat Cushion 		£
MAGWOSEATCUSH	No seat cushion	0
MAGCMRC	MPS Removable Seat Cushion Choose Between Seat Width 405 / 455 / 510 / 560 / 610. Choose Between Seat Depth 405 / 455 / 510. Choose also Upholstery in Section B - UPH For Rehab Seat Cushion	235
MAGPADDEDSC	Flat padded seat cushion Choose Seat Width And Seat Depth in Section B - Rehab Seat Sizes. Choose also Upholstery in Section B - UPH For Rehab Seat Cushion	0

UPH For Seat Cushion		£
MAGSUFABRIC	Fabric	0
MAGSUVINYL	Vinyl	0

Backrest 		£
MAGREHABCANES	Rehab canes, choose MPS or rehab insert back later	0
MAGREHABFOLD	Rehab canes with fold forward, choose MPS or rehab insert back later	375

Backrest Recline Option		£
-------------------------	--	---

Recline Option for Rehab seat and Rehab canes		
MAG3920-S-3921-F	Fixed (90° - 105°), (90°, 97,5°, 105°)	0
MAG3924-S-3925-F	Manual Recline (90° - 170°)	395
MAG3928-S-3929-F	Power Recline (90° - 160°) Backrest Height Minus Backrest Width Must Exceed 50 mm	1175
MAG3932-S-3933-F	Anti-shear power recline (95° -155°) Backrest Height Minus Backrest Width Must Exceed 50 mm, not with SW 305, 330, 355	1875

Recline option for Rehab seat and Rehab canes, Foldable to the front		
MAG3922-S-3923-F	Fixed (90° - 105°), (90°, 97,5°, 105°)	0
MAG3926-S-3927-F	Manual Recline (90° - 170°)	395
MAG3930-S-3931-F	Power Recline (90° - 160°) Backrest Height Minus Backrest Width Must Exceed 50 mm	1175
MAG3934-S-3935-F	Anti-shear power recline (95° -155°) Backrest Height Minus Backrest Width Must Exceed 50 mm, not with SW 305, 330, 355	1875

Insert Backrest		£
MAGBACKRESTINSERT	Magic Rehab insert backrest Only Available With MAGREHABCANES And MAGREHABFOLD	535
MAGMPSINSERT	MPS insert backrest Only with MTO fabric upholstery (except Standard), Choose Size in MPS Backrest - Sizes Section	0

Please Fill In Only If MPS Insert Backrest Chosen

MPS Insert Back - UPH		£
MAGFABRICBR	Fabric	0
MAGVINYLBR	Vinyl	0

MPS Insert Backrest – Sizes		£
-----------------------------	--	---

Backrest width Backrest Height Minus Backrest Width Must Exceed 50 mm		
MAGBW405MM	405 mm	0
MAGBW455MM	455 mm	0
MAGBW510MM	510 mm	0
MAGBW560MM	560 mm	0
MAGBW610MM	610 mm	0
Backrest height Backrest Height Minus Backrest Width Must Exceed 50 mm		
MAGBH535MM	535 mm	0
MAGBH585MM	585 mm	0

Please Fill In Only If Rehab Insert Backrest Chosen

Rehab Insert Backrest Cushion		£
MAGBRMRJ	Junior (350 mm height)	0
MAGBRMRJXS	Junior (Xtra tall, 410 mm height)	0
MAGBRMRMS	Medium (Short, 430 mm height)	0
MAGBRMRM	Medium (480 mm height)	0
MAGBRMRMXT	Medium (Xtra tall, 530 mm height)	0
MAGBRMRLS	Large (Short, 430 mm height)	0
MAGBRMRL	Large (480 mm height)	0
MAGBRMRLXT	Large (Xtra tall, 530 mm height)	0

Rehab Back Insert UPH		£
MAGFABRICBR	Fabric	0
MAGVINYLBR	Vinyl	0

Rehab Back Insert Bolster Size		£
For Magic Rehab back cushion large/medium		
MAGBOLSTERS	Shallow (110 mm depth)	0
MAGBOLSTERSTD	Standard (160 mm depth)	0
MAGBOLSTERD	Deep (220 mm depth)	0
MAGNOBOLSTER	No bolster	0

Rehab Back Insert Lumbar		£
MAGLUMBARS	Shallow (12 mm insert)	0
MAGLUMBARSTD	Standard (25 mm insert)	0
MAGLUMBARD	Deep (35 mm insert)	0
MAGNOLUMBAR	No lumbar	0

Insert Back Clamps		£
MAG3495	Standard, fixed clamps	0
MAG3492-PI	Quick Release	0

Push Handles		£
MAG0057-PI	MPS Push Rail	135
MAGPRHA	MPS Height adjustable handles	135
MAG0056-WOKOR-2	MPS Peg handles	135
MAG46354636	Rehab peg handles	135

Backrest Bag Hooks		£
MAG5724	MPS Bag hooks (pair) on backrest	35
MAG5723	Rehab Bag hooks (J-shaped) on spreader bar	35

Headrest		£							
MPS Headrest									
MAG2770VI	MPS Headrest, MPS Mount, Fabric to match backrest Mounting As For Car Seat Headrests, With Fabric Cover	270							
Rehab Headrest									
MAGREHABSTD	Magic Rehab Standard	160							
MAGREHABWIDE	Magic Rehab Wide	220							
Headrest Mount (for REHAB headrests only)									
MAG4113	2 AXIS (single stem, 2 pivot points)	165							
MAG4114	3 AXIS (single stem, 3 pivot points)	190							
MAG4115	MULTI AXIS (single stem)	215							
Rehab Headrest		£							
MAGFABRICHR	Fabric	0							
MAGVINYLHR	Vinyl	0							
Armrest – Types		£							
MAG3572	Pair - standard armrest, single stem	0							
MAG3994	Pair - standard armrest, double stem	160							
MAGFLEXIARM	Pair - flexi armrest	375							
MAG6503-6502	Pair - flexi armrest with locked position	550							
MAG3331-6500	Pair - flexi armrest with locking post Locking Post Attached To Seat Rail	600							
MAG3727-6501	Pair - flexi armrest with locking post and drop down	700							
Armrest – Heights									
MAGAH220MM	220 mm	MAGAH260MM	260 mm	MAGAH300MM	300 mm	MAGAH340MM	340 mm	MAGAH380MM	380 mm
MAGAH230MM	230 mm	MAGAH270MM	270 mm	MAGAH310MM	310 mm	MAGAH350MM	350 mm	MAGAH390MM	390 mm
MAGAH240MM	240 mm	MAGAH280MM	280 mm	MAGAH320MM	320 mm	MAGAH360MM	360 mm		
MAGAH250MM	250 mm	MAGAH290MM	290 mm	MAGAH330MM	330 mm	MAGAH370MM	370 mm		
Armpad – Style		£							
MAGAP75-WOKOR-2	75 mm armpads	0							
MAG3971-3972	MPS armpads	55							
MAG3712	MPS flip-out armpads	80							
MAGAP50-2	50 mm armpads	160							
MAG6711	Waterfall armpads	160							
MAG3775	Gutter armpads	160							
MAG6750	Gutter armpads with elbow support pad	160							
MAG6751	Gutter armpads with hand pad	160							
MAG6752	Gutter armpads with hand pad and elbow support pad	320							
Armpad – UPH		£							
MAGAPVINYL	Vinyl	0							
MAGAPFABRIC	Fabric	0							
Armrest Bag Hooks		£							
MAGHOOKA2RH	Right Side Armrest baghooks	35							
MAGHOOKA2LH	Left Side Armrest baghooks	35							

Footplate Height

MAGLL300MM	300 mm	MAGLL340MM	340 mm	MAGLL380MM Not with MAG6552	380 mm	MAGLL420MM Not with MAG6552	420 mm
MAGLL310MM	310 mm	MAGLL350MM	350 mm	MAGLL390MM Not with MAG6552	390 mm	MAGLL430MM Not with MAG6552	430 mm
MAGLL320MM	320 mm	MAGLL360MM	360 mm	MAGLL400MM Not with MAG6552	400 mm		
MAGLL330MM	330 mm	MAGLL370MM	370 mm	MAGLL410MM Not with MAG6552	410 mm		

Legrest – Choice

		£
MAGCPLEGREST	Centre Post legrest	0
MAGCPOWER	Centre Mount Power elevating legrest Not Available With Fixed Tilt MAG3502-KIT	1500
MAG3528	Swingaway legrest	0
MAG3527-ME	Manual Elevating Swinagaway legrest	425
MAG3526-ME	Power Elevating Swingaway legrest Not Available With Fixed Tilt MAG3502-KIT	850

Footplate

		£
MAG4055FIX	Centre Post: 1 piece fixed footplate (footplate width 325 mm)	0
MAG4055	Centre Post: 1 piece flip up footplate (footplate width 260 mm)	0
MAG4055-2	Centre Post: 2 piece flip up footplate (footplate width 260 mm)	0
MAG4755	CMPE: 1 piece flip up	0
MAG4755-2	CMPE: 2 piece flip up	0
MAGFLIPUP	2 piece flip up Only For Swing Away Legrests	0

Footplate Options

MAG4312	Swingaway footplate lock Only For Swing Away Legrests	55
MAG4260	Abduction block (neoperene)	55

Calf Support

		£
MAGCALFSUPPORT	Calf support standard Selection not necessary with MAGCPOWER, Is always included with MAGCPOWER	135

Knee Pads

		£
MAG4247	Knee pads	120
MAGKPCMPE	Knee pads for CMPE	120
MAG4660	Knee pads for non-elevating swingaway legrests	120
MAG4104	Knee pads for elevating swingaway legrests	120

Lap Belt

							£
MAG4973-MPS-4974-R	Standard Posture belt	500 mm	600 mm	700 mm	800 mm		0
MAG4973-4974-PBPAD	Standard Posture belt with 2 posture belt pads	500 mm	600 mm	700 mm	800 mm		100
MAG6540	4 point seat belt (PB4P)	500 mm	600 mm	700 mm	800 mm		180
MAG6540-2XPBPAD-4	4 point seat belt with 2 posture belt pads	500 mm	600 mm	700 mm	800 mm		1850
MAG4730	Retractable lap belt (for rehab seat)						110
	MAGLBLH					Left sided buckle	0
	MAGLBRH					Right sided buckle	0
MAG4731	Retractable lap belt (for MPS seat)						110
	MAGLBLH					Left sided buckle	0
	MAGLBRH					Right sided buckle	0

Thoracic Support (pad and mount)		£
MAG1506	Adjustable brackets (MPS back)	530
MAG4943	Adjustable brackets (MPS insert back)	530
Hip Pads - In Vinyl		£
MAGHPKIT	Hip pads on adjustable brackets (pair)	270
Harness		£
MAGHARP	Harness pediatric	160
MAGHARQ	Harness standard	160
Controller – R-Net		£
MAGPG-JSMLEDL	Basic LED JSM (no colour screen)	0
MAGPG-JSM-PI	Colour Screen JSM	160
MAGPG-JSM2IRL	Colour screen JSM2 + IR only	315
MAGPG-JSMIRBT	Colour screen JSM2 + Bluetooth + IR	530
Display Language Selection		£
MAG119001	English	0
Control+5 Module		£
MAG158709	CTRL+5 Switch Module (button)	320
MAG158708	CTRL+5 Switch Module (Toggle style)	320
Control+5 Location		£
MAG5430	On the armrest, same side as joystick	0
MAG5431	On the armrest, opposite side to joystick	0
MAG5432	Attendant (rear mounted)	100
Joystick Knob Option		£
MAGKNOBSTD	Standard joystick knob	0
MAGPG-KS80	Stick S80	80
MAGPG-KCHIN	Chin C38	80
MAGPG-KB32	Ball B32	80
MAGPG-KB38	Ball B38	80
MAGPG-KM48	Mushroom M48	80
MAGPG-KTBAR	T-Bar T80	80
MAGPG-KSB70	Soft Ball SB70	80
Joystick Position		£
MAGJSRIGHT	Joystick Right	0
MAGJSLEFT	Joystick Left	0
MAGPGCYPGC2.03725	Interchangeable to Left and Right	215
MAG3329JSM-RH	Mount JSM as an attendant control - Right (Rehab back)	315
MAG3329JSM-LH	Mount JSM as an attendant control - Left (Rehab back)	315

Swingback Options		£
MAG110050	Joystick Mounting, Fixed	0
MAG3083-KIT-OS	Swingback joystick to the outside (with magnets)	0
MAG3083-KIT-IS	Swingback joystick to the inside (with magnets)	0
MAG3082-KIT-OS	Swingback joystick to the outside (with locking latch)	215
MAG3082-KIT-IS	Swingback joystick to the inside (with locking latch)	215

Joystick Mounting Position		£
MAG4980	Height Adjustable Joystick Stem - for standard joystick stem	105
MAG3713	Height Adjustable Joystick Stem - for swingback joystick stem	105
MAG3729STD	Quick release stem (for standard joystick stem only)	0
MAG3729SB	Quick release stem (with swingback)	0

Speciality Controls		£
MAG5436	Omni2	1170
MAGPG-OMNISCANNER	Omni2 Enhanced 6-way scanner (single switch feature)	1170
	Omni2 location	
	MAGOMN1LH: Omni2 location, Left	MAGOMN1RH: Omni2 location, Right

Attendant Control		£
MAG3330AM-KIT	Attendant control module for MPS back	290
MAG3329AM-KIT	Attendant control module for Rehab back	290
	Attendant Control	
	MAGACLH: Attendant Control, Left	MAGACRH: Attendant Control, Right

Accessories		£
MAG4471	RAM® Cup holder	55
MAG4627	Magic Mobility Steel cup holder Drink Holder (Always Mounted On The Opposite Side Of The Control Panel/Table)	105
MAG4478	Magic Mobility Steel cup and phone holder (combined)	150
MAG4705	Accessory charger (1x12V + 2xUSB) - side rail mounted, same side as joystick	160
MAG6629	Accessory charger (2xUSB) - side rail mounted, opposite to joystick	160
MAG4801	Accessory charger (2xUSB) - armrest mounted, same side as joystick	160
MAG6594	Accessory charger (2xUSB) - armrest mounted, opposite to joystick	160
MAG5360	Loud horn	55
MAG4834	A4 Accessory bag (mounted on the backrest)	105
MAG4834RH	A4 Armrest bag (mounted on the right armrest)	105
MAG4834LH	A4 Armrest bag (mounted on the left armrest)	105
MAG4834-2	A4 Armrest bags (mounted on the left and right armrest)	210
MAG5320RHI	A5 Armrest bag (mounted internally, right armrest)	105
MAG5320LHI	A5 Armrest bag (mounted internally, left armrest)	105
MAG5320RHE	A5 Armrest bag (mounted externally, right armrest)	105
MAG5320LHE	A5 Armrest bag (mounted externally, left armrest)	105
MAG48345320I	A4 & A5 Armrest bags (A4 joystick side, A5 internal facing)	210
MAG48345320	A4 & A5 Armrest bags (A4 joystick side, A5 external facing)	210
MAG5699	Additional Armrest Receiver needed if more accessories should be mounted	35
MAGCOVER	Wheelchair cover	95

Clear Tables		£
MAG2677	Swingaway small table, 320 x 240mm not with MAGPGCYPGC2.03725, interchangeable joystick mount	335
MAG7700	Swingover table with magnet	335


Dimension Matrix MPS-SEAT

◀ BACK



Dimension Matrix REHAB-SEAT

◀ BACK



Seat Cushion For Rehab Seat[◀ BACK](#)

MAGCMRC

Backrest for Rehab Seat Pan[◀ BACK](#)

MAGMPSBACKREHAB

MAGIC REHAB INSERT SIZE




Chest Width 'C'	Height 'Y' (mm)
330-470mm	
Large (Xtra tall)	530
Large	480
Large (Short)	430
280-420mm	
Medium (Xtra tall)	530
Medium	480
Medium (Short)	430
260-400mm	
Junior (Xtra tall)	410
Junior	350

- Bolsters can be mounted upside down for increased armpit support
- Bolster heights can be offset up and down, up to 25mm between each side (asymmetric)

© Magic Mobility 2016

REHAB SEAT AND BACK SIZING



© Magic Mobility 2016

Colors

[◀ BACK](#)

Drive Wheels

[◀ BACK](#)

XT4



BLACK PNEUMATIC

OFF-ROAD

- Soft ride with good shock absorption
 - Best traction on all surfaces
 - 700mm base width
 - Tyre Liner available

Special Tyre Options

◀ BACK

MAGIC 360 PUNCTURE RESISTANCE

**NEVERFLAT**

- Offers good puncture protection for pneumatic tyres


Liner fits in between tube and tyre, providing a robust protective barrier

**PUNCTURE RESISTANT KIT**
[WITH PNEUMATIC OFF-ROAD & CROSSOVER TYRES]

- Offers comprehensive protection against glass, rocks, thorns, etc.
- Produced with aramid fibre and natural rubber
- Filled with Neverflat to give extra protection
- Can be retrofitted to all Magic Mobility chairs with knobby tyres

Headrest

◀ BACK

**MPS HEADREST**

ARMREST TYPE (1)

**STANDARD**

Note: 25mm of extra adjustment is added to bottom


**DOUBLE STEM**
For bariatric users**FLEXI**

Flip-up, paralleling


Note: Depending on requested armrest height the Flexi armrest will be assembled using slightly different parts. Therefore the range for each armrest will not be the full 220-440mm

Requested Height	Actual Range (Flexi on MPS back)	Requested Height	Actual Range (Flexi on Rehab back)
220-330mm	220-370mm	220-330mm	180-360mm
340-370mm	250-400mm	340-440mm	275-440mm
380-440mm	330-440mm		

ARMPAD ADJUSTMENTS (ALL ARMREST TYPES)







FLEXI ARMREST ADJUSTMENTS



ARMPADS



Allows the arm pads to be in tight to the user while still allowing flip up (Flexi arms only)




If MPS, 50mm, 75mm, gutter or waterfalls armpads are selected, the length of armpad fitted will generally be based on the seat depth. See below table for details. If specific armpad length is required please state this on order form

SEAT DEPTH ▶	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
50mm Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	← Armpad length
75mm Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
Gutter Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
Waterfall Arm Pad	255	295	305	305	305	355	355	355	380	380	380	380	380	
MPS Arm Pad	-	-	-	-	355	355	405	405	460	460	480	480	480	← Armpad length

Joystick Mounting

◀ BACK



Swingback Joystick Release Options

Magnet	Latch
- Requires more force but no dexterity to release - Can potentially release if wheelchair is being driven erratically or over rough terrain - Releases more easily if impacted	- Requires minimal force but some dexterity to release - Will stay locked at all times unless latch is pushed - Less likely to release if impacted, potentially causing more damage - Can be locked in fully swing-back position (optional) - Custom levers can be screwed to latch to make releasing it easier for clients with less hand movement



JOYSTICK STEM OPTIONS



65mm downwards adjustment from default position. Adjust height by releasing star knob



65mm downwards adjustment from default position. Adjust height with 3mm Allen key



HEIGHT ADJUSTABLE JOYSTICK STEM



QUICK RELEASE MOUNT

Joystick Knobs

◀ BACK



FOOTPLATE SELECTION



1 PIECE FLIP UP



2 PIECE FLIP UP





1 PIECE FIXED

Fastening bolts can be removed to allow footplate to flip up if required in the future



2 PIECE FIXED

Can adjust each footplate to a different angle




SHORT FOOTPLATES




© Magis Mobility 2016

Legrest MK2 Footplate Lock




Knee Pads


◀ BACK



**Centre mount legrest
(fixed and power elevating)**
3 axes of adjustment




Power or gas swingaway legrest
Fixed kneepad attached to legrest pivot



Standard swingaway legrest
3 axes of adjustment

Rear Boot Plate

◀ BACK



360

Lift up rear and pull
back to remove

Cup Holder

◀ BACK

Gradual configuration allows cup to stay flat when
tilting and when driving on bumpy surfaces



Internal diameter at
bottom is 84mm

If swingaway legrests are fitted the cup holder
mount bracket may be offset higher to avoid top of
legrest (depending on armrest height)



Note: Previous cup holder can be fitted if required.
See option on order form.



Armrest Bag

◀ BACK



Mounts under armrest
(unless otherwise requested)

If swingback joystick is selected it is generally more suitable to fit accessory bag to opposite side

© Magic Mob

Table Options

◀ BACK

TABLE OPTIONS

Clear polycarbonate table shown below. Black HD plastic is also available (typically without folded edges)

Unless a drawing has been supplied with specifications for cut out sizes for stomach and joystick, the slide in and swingover tables will be provided as a kit

Table stem slides into opposite side to joystick




**SLIDE IN
LARGE TABLE**

Table rests on armpad opposite to stem



**DOUBLE STEM
LARGE TABLE**

For chairs with attendant control ONLY (no joystick)



SWINGAWAY SMALL TABLE
Width 320 x Depth 240



SWINGOVER LARGE TABLE
Table rests on both arm pads

© Magic Mobility 2016

Comments

--

XT4

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day

Month

Year



Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel: +49 (0) 7253-9800

Fax: +49 (0) 7253-980111

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.de