

Otto Bock®

QUALITY FOR LIFE



Foot Caps

Enhanced stability and safety

The new Foot Caps fixation sandals from Otto Bock assure enhanced stability and safety. They are easy to install and suitable for all Otto Bock wheelchairs and mobility bases. Length-adjustable hook-and-loop closures facilitate optimum adaptation to the desired size. The anti-slip coating offers excellent support and stability while the protective back keeps the foot in the right position throughout the day.

Otto Bock HealthCare Plc

32 Parsonage Road, Englefield Green, Egham, Surrey, TW20 0LD - Tel: 01784 744900 - Fax: 01784 744901 - Email: bockuk@ottobock.com - Website: www.ottobock.co.uk

Foot Caps

Enhanced stability and safety




Technical data


	S (476B52=SK010)	M (476B52=SK015)	L (476B52=SK020)
Dimensions (LxW in mm)	190 x 89	230 x 114	255 x 140
Material	PA 6 30 % GF	PA 6 30 % GF	PA 6 30 % GF
Sole thickness	10 mm	10 mm	10 mm
Buckle width	32 mm	32 mm	32 mm
Rear edge height (heel)	35 mm	40 mm	45 mm
Weight	0.20 kg	0.28 kg	0.38 kg
Allowable body weight	75 kg	125 kg	180 kg
Available colours	Black	Black	Black



Size S



Size M



Size L

Otto Bock[®]

QUALITY FOR LIFE