

Invacare® | **Propad® Original Cushion**



**Designed for comfort, support and pressure reduction**

The *Invacare Propad Original* cushion has been designed to give comfort and pressure reduction. The *Propad Original* cushion is suitable for those considered 'At Risk'\* of developing pressure ulcers.



**Yes, you can.®**



**Features and Options**

**Deep profiled foam cushion, designed to deliver pressure care**

The *Invacare Propad Original* Cushion is made from a layer of castellated foam. Each castellation moves independently to deliver enhanced support and pressure area care to the client, as well as good air circulation, reducing the potential for sweating and associated skin maceration.


**Contoured seating surface**

Offers improved stability and postural support; castellations reduce shear and friction forces on tissue.



**Cover helps prevent a build up of shear and friction forces against the client**

Cover is multi stretch, water-resistant and vapour-permeable. Cover features an anti-slip base.



**Available in a wide range of sizes**

To view the cushion size range, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

**Technical data**



Wide range of sizes



Height: 102 mm / 4"



Weight: 1 kg



Max User Weight: 110kg / 17.3 stone



Cleaning: 80°

**Propad Original Cushion**

Manufactured to European Standards:  
Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5)  
Cover: BS 7175 (Crib 5)

**Warranty 1 Year**

INVACARE® Limited  
Pencoed Technology Park - Pencoed - Bridgend CF35 5HZ  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)

INVACARE® Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County Dublin  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

**Colours**



Blue



Black Towelling

© 2010 Invacare International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Propad Original Cushion  
UK - 10/2010 - 1505219

\* Essential nursing care is pivotal in pressure ulcer prevention. This product will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.



**Yes, you can.®**