



Side Rail Height Extension



Bed Extension



Mattress Extension



Bariatric Grab Handle

Nursing Care Beds / Bradshaw Bariatric Low Nursing Care Bed

MATTRESSES

See pages 59-89

Bradshaw Bariatric Low Nursing Care Bed:

1275/BAR/LOW/LOAK/S

Key Features:

Robust and sturdy electrically profiling Bariatric bed for use in the Community, Nursing and Residential Care markets.

Quick and easy to assemble, even within limited spaces.

Low platform height range of 25.9cm to 62.5cm.

Lockable handset controls all electrically operated functions.

Four braked castors for added safety.

Electrically operated rising backrest and kneebreak functions.

Infinitely variable tilt and reverse tilt positioning.

Attractive domestic appearance with high quality wooden head and foot ends and full length integral wooden side rails.

There is no provision for use of an integrated lifting pole; instead, for bariatric clients, we would recommend the use of an appropriate hoist.

The four section bed can be dismantled and stored on a transport stand.

3 year warranty (1st year parts & labour, 2nd and 3rd year parts only).

Significant increase in engineering specification means that the bed can support a maximum user weight of up to 320kg (50st).

Accessories:

Elland Bradshaw Bariatric Grab Handle	1215/GRAB
Bradshaw Bariatric Bed Extension Kit (Includes longer side rails)	1275/BAR/EX/LOAK
Cromer Cantilever (Bed Cradle)	1369
Spare Side Rails (Bradshaw Bariatric)	1275/BAR/SR/LOAK
Spare Side Rail (Single)	1275/SIDE/LOAK
Bradshaw Bariatric Side Rail Height Extension (Pair)	1275/5/BAR/LOAK
Bradshaw Side Rail Pads (Full Length)	1316/PS/74/CREAM
Bradshaw Bariatric Spare Handset	CONTROL/046
Bradshaw Bariatric High Side Rail Pads (Pair)	1316/PS/75/HI/CREAM

Height (cm)	Width (cm)	Length (cm)
25.9 to 62.5 (to mattress platform)	136 (incl. side rails)	225
118 (overall bed height at highest point)		
Product Weight (kg)	147	

For suitable surfaces (mattresses) please see the Selector on pages 86-89.
For easy selection side rail and bed accessory tables, please see pages 41 and 44-45, respectively.