

SmartDrive MX2+

Connect and go



permobil

Energy to do more

Pushing a manual wheelchair can put stress on the upper body and lead to serious issues in the future. SmartDrive gives you freedom to do more. A lightweight, rear-mounted power assist device that moves seamlessly with your chair, even when not in use. **Simply connect and go.**

- Compact and convenient for easy attachment and removal within the existing footprint of the chair
- Start, stop and change speed on-the-go
- Available with multiple control options including wearable and wired options so you can power your push, your way
- Easy to transport
- Range of up to 19.8km



Simple

SmartDrive provides manual wheelchair users more control over their experience than ever before. Discover our range of control options to better support your individual needs.



	SpeedControl Dial	SwitchControl	SwitchControl mono jack	SmartWatch*
Provided by Permobil	✓	✓	✓	
Wired option	✓	✓	✓	
Plug & play	✓	✓	✓	
Set up within 5 min	✓	✓	✓	
Programmable setting	✓	✓	✓	✓
No charging required	✓	✓	✓	
Works with SpeedControl	✓			✓
Works with SwitchControl		✓		✓
Works with SwitchControl mono jack			✓	✓
No Wi-Fi required	✓	✓		
Increase or decrease speeds without stopping	✓			

*Available on Apple Watch and Samsung Galaxy

Specifications

Dimensions	Length: 389 mm Width: 141 mm Height: 242 mm
Product weight	5.7 kg Drive Unit
User weight	14 to 150 kg
Driving range*	Up to 19.8 km Level/No Resistance
Motor power	250 W Brushless DC
Maximum speed	6 km/h
Operating temp	-25°C to 50°C
Chair type compatibility	Rigid, Folding, Tilt-In Space, One-Arm Drive and others
Drive wheel diameter	501, 540, 559, 590 mm
Warranty	2-Year Limited Warranty
Bluetooth® LE	FCC QOQBT113 ID 5123A-BGTBLE113 IC ID

All figures shown are for illustration purposes only and actual products may vary. Wear OS by Google, iPhone, Android, Bluetooth and Google Play are trademarks or registered trademarks of their respective owners.

