



# PUMA 20

SIMPLY GOOD DRIVING



# PUMA20



# SIMPLY GOOD DRIVING

You expect your powered wheelchair to take you comfortably wherever you want to go. And offer you all the essential functions you need for an active daily life. A modern mobility product that fits your lifestyle. This is fully reflected in the Puma 20. Robust, practical and affordable, this power chair offers simply good driving in every situation - and many optional extras.

**// Best driving performance in its class:**

Enjoy maximum independence indoors and out thanks to a technically smart and robust design, good suspension and powerful motors.

**// Adaptable to suit your needs:**

Combine the Sedeo Lite seat frame with your choice of Sedeo comfort seating options and electrical functions for enhanced comfort and safety.

**// Distinguishing design for a pleasant and carefree ride:**

Highly functional and durable with an eye for aesthetics - enjoy comfortable and safe driving in style.

**// Front-wheel or rear-wheel drive:**

Choose the drive base that fits your driving preferences. And swap around if required.





# MAXIMUM INDEPENDENCE FOR AN ACTIVE LIFE

Designed for intensive use, the Puma 20 delivers the best driving performance in its class. Naturally, it's easy to drive this powered wheelchair with your choice of renowned electronics. Enjoy a smooth and stable ride wherever you go thanks to the robust base, good suspension – with rubber shock absorbers - and comfort seating. More freedom outdoors? Choose a higher top speed and larger range with the 74 Ah battery option. High manoeuvrability coupled with a compact design and low seat height make for great use in and around the house.

**// Front-wheel drive or rear-wheel drive base**

Choose front-wheel drive for even more indoor manoeuvrability. Rear-wheel drive ensures maximum outdoor performance. Changing needs? Your Puma 20 power chair can be easily converted from front-wheel drive to rear-wheel drive or vice versa.



# ALL THE ESSENTIALS - AND MORE

With the Puma 20 you get simply good driving - with the option of much more... Truly modular in design, the Puma 20 offers you the flexibility to easily 'upgrade' your chair to suit your personal needs. Combine the robust Sedeo Lite seat frame with your choice of high-quality Sedeo seat cushion, backrest, armrests and legrests. And choose electrical 25° seat tilt and lights for enhanced seating comfort and safety.





DENHAM | MAISON SCOTCH  
DEPARTMENT | KING | FOL  
A | CARHARTT | DEPT | M  
IRAAT | MUSSY | 8 MILLIMETER  
DIESEL | DENHAM | IC. REG. | LE

# DISTINGUISHING DESIGN FOR PLEASANT AND CAREFREE DRIVING

Sunrise Medical has a reputation for quality and safety. Do not expect anything less from your Puma 20. Our uncompromising approach to performance, comfort and ease of use has resulted in a technically smart, modular design with lots of possibilities. Thanks to the robust base and good suspension, this chair offers a surprisingly good driving experience. While creating a highly functional and durable wheelchair, however, we've certainly not overlooked the importance of aesthetics. Enjoy comfortable and safe driving in style!



# TECHNICAL DATA

## PUMA 20 // POWERED WHEELCHAIR

<b>Max. user weight:</b>	136 kg	<b>Max. reversing width:</b>	1200 mm (FWD); 1700 mm (RWD)
<b>Seat width:</b>	400 – 510 mm	<b>Max. safe slope:</b>	6° (FWD); 8° (RWD)
<b>Seat depth:</b>	420 – 500 mm	<b>Max. kerb climbing:</b>	60 mm (FWD); 80 mm (RWD)
<b>Seat height:</b>	from 405 mm	<b>Wheel options:</b>	12.5" drive wheels 8 or 9" castors; 14" drive wheels & 10" castors
<b>Seat tilt:</b>	0 - 25°	<b>Colour choice:</b>	Anthracite grey / Platinum white base
<b>Backrest height:</b>	480 mm	<b>Max. range (ISO 7176-4)</b>	42.5 km (with 74 Ah batteries & 6 km/h motors, dependent on distance profile, user weight and electrical options.
<b>Backrest angle:</b>	92 / 97 / 102 / 107°	<b>Wheelchair weight:</b>	from 106 kg (including 40 Ah batteries)
<b>Total width:</b>	610 mm (12.5" drive wheels); 655 mm (14" drive wheels)	<b>Crash tested:</b>	Yes (approved for transport in a vehicle according to ISO 7176-19)
<b>Total length:</b>	1160 mm (FWD); 1095 mm (RWD)	<b>Transport without user:</b>	Fix the wheelchair safe in the vehicle (e.g. car or plane)
<b>Total height:</b>	from 1030 mm	<b>Intended use and environment:</b>	For people with limited mobility; indoor and outdoor use. Wheelchair type class B
<b>Max. speed:</b>	6 / 10 km/h		
<b>Battery size:</b>	40 / 50 / 60 / 74 Ah		
<b>Electronics:</b>	VR2 / R-net controller		

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.eu](http://www.SunriseMedical.eu) or alternatively is available on request in large text.



### Sunrise Medical HCM B.V.

Vossenbeemd 104  
5705 CL Helmond  
The Netherlands  
T: +31 (0)492 593 888  
F: +31 (0)492 537 931  
E: [customerservice@sunrisemedical.nl](mailto:customerservice@sunrisemedical.nl)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)