

Uno M Mobile Lift

For Use in Home
and Care Facilities



FEATURES AND BENEFITS¹

- **Large lifting interval** accommodates a variety of lifting situations
- **FlexLink** facilitates a smooth lift motion and a greater clearance to the lift mast
- **Manual and electric emergency lowering mechanisms** are designed to help support operation for caregivers
- Compact design offers a **range of access for lifting** throughout the home and care facilities
- **Intuitive user display** includes lift count, battery usage and overload indicator
- **Optional quick release hook system** helps caregivers change a slingbar or install a scale
- **Optional Li-ion battery**
- **Extensive portfolio of accessories** includes slings, slingbars and scales



FlexLink facilitates a smooth lift motion and a greater clearance to the lift mast



Manual emergency lowering



Controls are standardised across Baxter mobile lifts

TECHNICAL SPECIFICATIONS

Maximum load	182 kg (400 lb)	
Weight	Uno M lift (excl sling bar, battery, charging cable) 36.9 kg (81.2 lb)	
Height	1333 to 2002 mm (52.5 to 78.8")	
Protection class (Lift & Pendant)	IP 24	
Slingbar options	Standard Universal SlingBar 450 (450 mm)	Optional Universal SlingBar 350, 600 Universal TwinBar 670 Universal SideBars 450 Sling Cross-bar 450, 670
Base	Height: 104 mm (4.7") Length: 1220 mm (48") Width: 681 - 1141 mm (26.8 - 44.9")	
Base-width expansion	Electric	
Casters:	Front 75 mm (3") twin casters	Rear 75 mm (3") twin casters with brakes
Material	Powder coated steel	
Emergency lowering	Manual and electric	
Lifting Speed (no load)	33 mm/s (1.3 in/s)	
Maximum lifting height	1783 mm (70.2")	
Lifting interval	1297 mm (51.1")	
Sound Level (max load)	51.7 dB	
Batteries	Lead-acid gel, valve-regulated battery - 24 V 2.9 Ah Prod. No. 2006106 Li-ion, Lithium Iron Phosphate (LiFePO4) 25.6 V 3.3 Ah Prod. No. 2006110	Weight Battery lead-acid gel: 2.8 kg (6.2 lb) Battery Li-ion: 1.5 kg (3.3 lb)
Charger	Built-in	
Product number	2050020	

For more information, please contact your local Baxter representative.

References:

1. Instructions for Use

Uno M mobile lift is a general-purpose lift with the intended use in; health care, intensive care, rehabilitation, and home healthcare environments. It is an excellent aid in daily transfers of both adults and children, for instance, lifting to and from wheel chair, bed, toilet, bathtub, and floor.

Medical device (2017/745): Class I. Manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå, Sweden. These products are regulated health care products which, where required by applicable regulations, bear a CE mark. Baxter recommends you carefully read the detailed instructions for safe and proper use included in the documentation accompanying medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

Hillrom.Baxter.eu/en-GB

Baxter Healthcare Ltd

Spur A, 2nd Floor, Reading International Business Park, Basingstoke Road, Reading, RG2 6DH

Baxter and Uno are trademarks of Baxter International Inc. or its subsidiaries.

UKI-CS123-260001 (v1.0) 04/2026