

  
**zippie**<sup>®</sup>  
Explore your world!



**Youngster** 



*Fit for Adventures*





**Youngster<sup>3</sup>**

The lightest folding wheelchair with growth adjustability!

**Zippie Youngster 3:**

The new lightweight Youngster 3 offers integrated seat width and seat depth growth. The new frame design gives your wheelchair a stylish and unique look.

On top of that the Youngster 3 offers new, innovative and lightweight options using the latest technology.







**PATENTED**

Patented 2cm seat width growth via cross brace for growth (without new parts).



New lightweight fork and frame for maximum adjustability, improved looks and better manoeuvrability.

**- Lightweight package:**  
Lightweight fork, lightweight seat upholstery, lightweight anti-tips, Safari wheel lock



4cm seat depth growth via seat tubes for growth (without new parts).

Explore your world! Have fun, enjoy playing with your friends and be active

### Options

The new Youngster 3 offers a wide range of options for more comfort, freedom, fun and safety.



New selection of back upholsteries with a wide range of different versions and improved adjustability.



Lightweight seat upholstery for easy handling and transport.

### Twenty eight frame colours!

Select a colour for each component.



Modern side guard design in combination with the new Safari wheel lock for improved safety.


### Spoke guards

There are 32 fantastic spoke guard designs you can choose from. See our spoke guard catalogue for the whole range.






# Technical specifications



<b>Frame colour :</b>	28 colours available	<b>Weight of heaviest part:</b>	4.4kg
<b>Upholstery colour:</b>	Black / Vented (black & silver) / Alcantara (grey & brown)	<b>Product weight:</b>	from 8.5kg
<b>Max user weight:</b>	85kg (13.5 stone)	<b>Front wheel diameter:</b>	10, 12, 15cm (4", 5", 6")
<b>Growth: Seat width:</b>	+2cm (0.7") without parts	<b>Rear wheel diameter:</b>	51, 55, 61cm (20", 22", 24")
<b>Seat depth:</b>	+4cm (1.6") without parts		

## Frame Colours



It is best to select colours using the metal colour plates at your dealers and to use the order form to finalise your desired colour scheme

### Sunrise Medical Limited

High Street,  
Wollaston,  
Stourbridge  
West Midlands  
DY8 4PS  
Tel.: +44 (0) 845 605 6688  
Fax: +44 (0) 845 605 6689  
www.SunriseMedical.com



PB\_Youngster3\_EN\_GB\_REV2

### Find out more today!

Visit our website at [www.SunriseMedical.com](http://www.SunriseMedical.com) for more information or to request a demonstration.