

Otto Bock®

QUALITY FOR LIFE



C2000

Get There Faster

This is now very easy with the C2000 – the power wheelchair for outdoor use. An innovative steering unit ensures safety at speeds of up to 15 km/h, even on hills and uneven surfaces. The springs of the entire chassis along with the proper seat from the Otto Bock seating range provide maximum comfort on stretches as long as 60 km.

Don't think twice about driving outdoors - integrated lighting and a serialized carrier are included in the standard equipment of the front- or rear-wheel drive C2000.

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone +49 5527 848 1590 · Fax +49 5527 848 1676 · reha-export@ottobock.de · www.ottobock.com

C2000

Outdoor Wheelchair

Features

- Innovative force steering
- Seat raising of 40 cm (optional)
- Electrical tilting and back angle adjustment (optional)
- Various footplates for different areas of application/weights
- Carrier included
- Lighting included
- enAble50 steering, Otto Bock special steering available



Technical Data

Speed	6,10, or 15km/h
Battery capacity	86 Ah
Range	60 km
Electronic control	enAble50
Load capacity	Standard 160 (352 lbs), 200 kg (440 lbs), Special design <260 kg (572 lbs)
Obstacle height the wheelchair is able to overcome	150 mm
Climbing ability	22 %
Width	680 mm
Length with footplates	122 mm
Turning radius	940 mm
Seat width	Standard 38 / 48 cm
Seat depth	Standard 38 / 50 cm
Seat tilt	30°, electrical, continuously adjustable
Back angle	Standard -9 / 1 / 11 / 21°, 0 / 10 / 20 / 30°, electrical, continuously adjustable
Back height	55 cm
Lower leg length	40–55 cm
Seat height	<51 cm
Armrest height	25–36 cm

*Technical modifications reserved

Otto Bock[®]

QUALITY FOR LIFE