

#### Invacare Matrx Flo-tech Plus

The *Invacare Matrx Flo-tech Plus* cushion incorporates a gel sac which is located towards the rear of the cushion. The positioning of the gel sac accommodates the ischial tuberosities, helping to significantly reduce pressure in this vulnerable area. This delivers superior levels of protection for those patients identified as being at High Risk of developing pressure ulcers.

## Deep contoured foam base with gel sac at rear of the cushion

The positioning of the gel sac affords increased immersion and protection to the Ischial Tuberosities for users who are at High Risk of developing pressure ulcers.

#### **Contoured seating surface**

Offers improved stability and postural support. Also delivers a greater degree of weight distribution compared to that of a flat surface. The castellations of the cushion allow slight movement which in

turn reduces shear and friction forces. The deeper moulded foam offers greater protection to those who are considered clinically obese, or users who would benefit from a higher sitting position.

#### **Features and Options**



#### **Ingress Resistant Cover**

Cover is two-way stretch, waterresistant and vapour-permeable. Cover features an anti-slip base with carry handle.



#### Available in 80 sizes

The Invacare Matrx Flo-tech Plus is available in a range of 80 sizes, from  $36 \times 33$  to  $60 \times 55$  cm  $(360 \times 330 \text{ to } 600 \times 550 \text{ mm})$ .



### Integral Sag Compensator option

Option of sag compensator to accommodate the sag in wheelchair seats, and ensure the cushion is correctly presented.



# Ultra-Fresh™ added to all *Flotech* foam for antimicrobial and odour protection

Antimicrobial treatment that has been added to all *Flo-tech* foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.