ORDER FORM



ECO ASSIST base price £ 2695.00

ORDER	
QUOTE	

Company/Customer ID:		P.O. Number:	Date:	
Contact person:		Other delivery address/Company:		
Street address:		Street address:		
City and Post Code:		City and Post Code:		
TEL:	E-Mail:	Remark:		



User Information:	Requested delivery date:
User name as reference:	S/N1 (Rehasense to fill):
Wheelchair Brand/Model:	S/N 2 (Rehasense to fill):
Front Tube Diameter:	S/N 3 (Rehasense to fill):
Distance between front tubes (center to center):	S/N 4 (Rehasense to fill):
Form of the front tubes: ROUND OVAL KÜSCHALL K2 MEYRA NEW/	OTHER



1. ECO ASSIST

Picture	Specifications	Price excl. VAT	Select the type	Qty
	 Motor: 350W Tyre 12½"x 2½" with city pattern Battery: Lithium-ion, 48V, 5.8Ah (278Wh) Docking: Manual only Wheelchair fitting width range: 210-480 mm 	REA12MMS00 £ 2695.00 ●	with standard handles	
	 Wheelchair & frame type compatible: rigid & folding frames/ fixed front or swing away Max user weight: 110 kg Max speed: 20 km/h Transport weight: 10 kg/ total weight: 15.45 kg 	REA12MMT00	with Tetra handles	

2. Arms for width adjustment

Picture	Picture	Description	Extra price excl. VAT	Opt.	Qty.
	第1414041411	Arm 100 mm with straight clamp holder (default - including ECO Assist) (for mounting widths 270-380 mm) Item number: EBC202A-BK (1 piece)	المحادية ما	\circ	
center to center	新有种的 (C.)	Arm 150 mm with straight clamp holder (also including ECO Assist) (for mounting widths 370-450 mm) Item number: EBC202B-BK (1 piece)	included		
	(14) (14) (14) (14) (14) (14) (14) (14)	Arm 180 mm with straight clamp holder (for mounting widths 430-480 mm) Item number: EBC202C-BK (1 piece)	£ 0.00		
		Arm 100 mm with 9.5 degrees angled clamp holder - RIGHT (for mounting widths 210/220-320/330 mm) Item number: EBC212A-BK (1 piece)	5000		
Noto	WWW001	Arm 100 mm with 9.5 degrees angled clamp holder - LEFT (for mounting widths 210/220-320/330 mm) Item number: EBC212B-BK (1 piece)	€ 0.00		

Note

The width must be measured from the center to the center of the front tubes as illustrated in the picture above. It should also be mounted and measured as high up on the front frame as possible and in the middle of a straight piece of at least 5-6 cm, where the clamps should be placed. Usually the measurement do not have to be 100% accurate, it is good enough with a +/- tolerance of 5-10 mm.

Any additional info / remark?

3. Clamps

Picture	Description	Extra price excl. VAT	Opt.
	Standard manual clamp, black – LEFT Item number: EBX035-BK		
	Standard manual clamp, black - RIGHT Item number: EBX036-BK	included	
	Manual clamp with safety lock, black - LEFT (new) Item number: coming soon	€ 0.00	
	Manual clamp with safety lock, black - RIGHT (new) Item number: coming soon	€ 0.00	
•	Manual clamp for OVAL tubes, black - LEFT (new) Item number: coming soon	6.0.00	
Coming soon	Manual clamp for OVAL tubes, black - RIGHT (new) Item number: coming soon	£ 0.00	
1 1	Manual clamp with 40 mm high claws, black - LEFT Item number: EBX03A-BK	5000	
	Manual clamp with 40 mm high claws, black - RIGHT Item number: EBX03B-BK	£ 0.00	
Any additional info / remark?			



4. Longer spindle for manual clamps

Picture	Description	Extra price excl. VAT	Optional	Qty.
	Longer spindle for manual clamps required for tube diameter over ø32 mm Item number: EBX058L (1 piece)	€ 0.00		
Important: Please note that you need 2 pcs per whe	eelchair)			

5. Clamp tubes

Picture	Description	Extra price excl. VAT	Optional	Qty.
	Standard clamp tube for ECO ASSIST, 100 mm Item number: EBC203 (1 piece)	included	\circ	
	Standard clamp tube for ECO ASSIST, 150 mm Item number: EBC213A-BK (1 piece)	£ 0.00	0	
Important: Please note that you need 2 pcs per who	eelchair)			

6. Speed settings

Cear	1	2	3	Reverse	
Customer-specific gear and speed details					
Predefined speed	Yes please use my default setting profile			\circ	
Note: Always drive within your ability to control the device to a safe stop. Using/driving the device at uncontrollable speeds is the user's responsibility and the manufacturer/dealer will deny/refuse any responsibility of accidents or damages to third party. Any additional info / remark?					

7. Tube protectors

Picture	Description	Extra price excl. VAT	Optional	Qty.
	Tube protector for Panthera wheelchairs Item number: EBCO2N (1 piece)	£ 30.00		
	Tube protector for Ø22 mm tube Item number: EBC06N (1 piece)	£ 30.00		
	Tube protector for ø25 mm tube Item number: EBC05N (1 piece)	£ 30.00		
Important: Please note that you need 2 pcs per wheelchair)				



8. Alternative throttle levers

Pictures	Description	Extra price excl. VAT	Optional
WALK ON OFF	Right-hand throttle lever – default Item number:	included	
ON OFF WALK ON OFF	Left-hand throttle lever Item number: EBE00M-LT	£ 0.00	
WALK ON OFF	Thumb throttle for PAWS - fits both sides Item number: EBE00M	£ 50.00	
Where you want it mounted?	Left side: Right side:		
	Standard tetra handles - horizontal Item number (set of 2 pcs): Item number: EBH00L Included as standard if REA12MMTOO is selected	£ 0.00	
Any additional info / remark?			

9. Throttle Grip Assist

Picture	Description	Price excl. VAT	Optional	Qty.
ON OFF ON OFF	Throttle Grip Assist for standard handle (right hand only) - BLACK Item number: EBH218-BK	£ 10,00		
	Throttle Grip Assist for standard handle (right hand only) - BLUE Item number: EBH218-BE	£ 10,00		
	Throttle Grip Assist for standard handle (right hand only) - RED Item number: EBH218-RD	€ 10,00		
	Throttle Grip Assist for standard handle (right hand only) - GREEN Item number: EBH218-GN	£ 10,00		
	Throttle Crip Assist for standard handle (right hand only) - YELLOW Item number: EBH218-YW	£ 10,00		



10. Handlebar Stem

Pictures	Description	Extra price excl. VAT	Optional	Qty.
	Standard with 74 mm fixed length Item number: EBH00G-A	included		
	Short - length adjustable 180-215 mm Item number: EBH00G	£ 50.00		
	Long - length adjustable 210-265 mm Item number: EBH01G	£ 50.00		
	Clamp for handlebar stem and 90 degree angle. May be necessary in specific situations. Item number: EBH200-A	£ 30.00		

11. Alternative tyre/pattern

Pic	tures	Description	Extra price excl. VAT	Optional
Default tyre pattern	Alternative tyre pattern	Available tyre for 12" PAWS 12 1/2" x 2.75" with off-road pattern, black Item number: 431210	£ 50.00	
Info: The standard tyre will be shipped to you together with your ECO ASSIST				

12. Want to upgrade the 5.8 Ah battery for the larger 11.6 Ah battery and pay the difference?

Pictures	Description	Extra price excl. VAT	Optional
ECO assist	High capacity Lithium-ion battery, 48V, 11.6 Ah, 556.8 Wh - 3.1 kg Item number: EBB0048T2A	£ 250,00	



13. Additional accessories

Pictures	Description	Extra price excl. VAT	Optional	Qty.
eccessist •	High capacity Lithium-ion battery, 48V, 11.6 Ah, 556.8 Wh – 3.1 kg Item number: EBB0048T2A	€ 399.00		
	Rear view mirror Ø100 mm - 1 piece, can be used on both sides Item number: EBH01B	£ 70.00		
	Smart tyre inflator Item number: EBE01D	£ 65.00		
	Phone holder 1 - black aluminium, 360 degree rotation, width range 60-95 mm, thickness of the phone max 11 mm Item number: EBH220	£ 35.00		
	Phone holder 2 - black aluminium, 360 degree rotation, width range 60-100 mm, thickness of the phone max 12 mm Item number: EBH221	£ 40.00		
Coming soon	Bag + holder and bracket for ECCO ASSIST Item number: coming soon	n/a		
	Set of weights for ECO ASSIST - 10 kg Item number: coming soon	n/a		
Coming soon	Footrest solution for ECO ASSIST. Designed for use with wheelchairs with removable leg supports. Item number: coming soon	n/a		
Coming soon	Transport bag for ECO ASSIST Item number: coming soon	n/a		

■ Always in stock ONot in stock at the moment. Ask for delivery time.

Note

Always drive within your ability to control the device to a safe stop. Using/driving the device at uncontrollable speeds is the user's responsibility and the manufacturer/dealer will deny/refuse any responsibility of accidents or damages to third party.

Disclaimer to be signed if custom speed has been set up We herewith advise, that using the PAWS with max. speed-settings higher than 4 Mph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the PAWS with speed-settings higher than 4 Mph may not be driven on public roads, bicycle-lanes and pavements.

You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the PAWS, speed settings higher than 4 Mph (8,12 and 18 Mph). When you're travelling with your PAWS, always check the specific legislation in the respective country.

When placing an order for a max, speed-setting higher than 4 Mph please confirm with your signature below, that you have understood the above information and that you only intend the usage of your PAWS on private property.

We herewith advise, that using the PAWS with max. speed-settings higher than 4 Mph is only permitted on public roads in certain countries and otherwise may only be used on private property.

please confirm with your signature below, that you have understood the above information and that you only intend the usage of your PAWS on private property herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.

Parts.	Contain of Contains
Date	Customer Signature