



## **Our Products**











Hi Law Seating		
<b>Hi-Low Seating</b> Hi-Low MPS	2	**
Multi-Positioning Seat Systems		
MPS	4	**
MPS Accessories	5	**
Customize your Size	6	**
MPS Car Seat	8	**
Soft-Touch® Sitters		
Sitter with Stationary Base	9	**
Sitter with Mobile Base	9	**
Floor Sitter	9	**
Sitter Only	10	*
Simple Seating Solutions		
Mobile Activity Tray	12	*
Height Right™ Chair	14	*
Height Right™ Accessories	15	**
Liners	16	*
Booster Seats	17	*
Portable Potty Seat	18	*
Soft-Touch® Therapy Products		
Therapy Wedges	19	**
Therapy Rolls	19	**
Push Chairs		
EIO	20	**
Jogger	22	**
Measurement Help	24	
Brand Story	25	

### **Positioning Support Levels**

- For children with minimal positioning needs: Child can sit on the floor independently, may have mild physical delays, displays good head and trunk control.
- \*\* For children with **minimal to mild** positioning needs:

  Child has some limited sitting balance, displays some independent control of head and trunk, has minimal restrictions of hip and knee range of motion, presents no significant issues with muscle tone (low or high).
- \*\*\* For children with **mild to moderate** positioning needs:
  Child requires support for sitting, muscle tone effects independent positioning (low or high), displays some control of head and trunk with support, child benefits from tilted position.

# **Special Tomato**®



#### **Soft-Touch® Material**

The Soft-Touch® material outer skin prevents the products from peeling or tearing and also prevents liquids from penetrating the core. It is cool and smooth reducing the possibility of children's skin sticking to the seating surface. The soft-to-the-touch inner core is a specially formulated foam that allows for children to sit comfortably for longer periods of time.



### **Safety Concern**

With over 35 years of manufacturing pediatric seating and positioning systems, our products are safe and easy to use while providing a latex-free contact surface. The contoured designs of our products provide a large surface area for distributing pressure which reduces the risk of skin shear and breakdown. Many of our products are FAA approved for airline use.



### **Design and Function**

Hand-crafted in upstate New York, the Special Tomato® products are designed with the needs of the user and care provider in mind. The superior Soft-Touch® materials are very durable and easy to clean. Our unique bent wooden bases deliver support and strength while the Soft-Touch® material allows our designers to subtlety integrate the positioning contours needed for a proper seating system. We stand behind the quality of our products by offering up to a 5 year warranty on most of our products.



### **Sales and Support**

Moorings Mediquip distribute Special Tomato in the UK and Ireland. Advice and assessments are available in all areas. Call 0800 031 6571 or email sales@mooringsmediquip.com.

1

# Hi-Low (MPS) Seating System

Height Adjustable Seating with Built-in Growth

NOW... More seat cushion choices from Early Intervention to Adult!

- Affordable Customized Seating
- Tilt-in-Space with Child in Seat
- Floor to Table Height Adjustment
- Contoured Cushion Support
- Tool Free Adjustments
- Ergo Harness Design
- Soft-Touch® Material Latex-Free



See page 7 for warranty details














- Ergonomic Push Handle
   Comfortable grip for all heights.
- 4 Tool Free
  Height adjustment of all cushions.
- MPS Tray 79060159 Dark Gray, Small 79060359 Dark Gray, Large Sold separately.

- 2 Wood Accents Added structure with a modern design.
- Footrest Angle, Depth and Height Adjustment Allows for optimal positioning.
- 8 Hi-Low Adjustment Raise and lower with child in seat.

- 3 **Tilt-in-Space**-30° to 10° range of tilt adjustment with child in seat.
- 6 Swivel Locking Casters
  All 4 easy lock 3" swivel casters
  for safety and control.
- Colorful Contoured Comfort Vibrant color choices and many cushion sizes available.



# Multi-Positioning Seat (MPS) 5 year warranty



#### Versatile Seating with Built-in Growth

NEW Cushion sizes available for both Small and Large MPS! Cushions grow the seating system from 9 months to young adult!

- Adjustable 5-point harness
- Hip-Flex® Seat Cushion
- Height adjustable Trunk and Head Cushions





**Adjustable Tray** Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.



Hip-Flex® adjustment provides a range of positions between 15° of hip flexion and 15° of hip

## Use the Special Tomato® MPS system everywhere!



At home, secures to most chairs by using attachment straps



On the school bus, transport the Special Tomato® in Push Chair (79200000, see next page)



At school, secure the Special Tomato® in a classroom chair, or use it with the Mobile Base as your all-in-one classroom chair



In an airplane, both sizes approved by FAA for air travel for seats purchased after April 1,2006





- 1 Mobile Floor Base 79070100 Small 79070300 Large
- Optional Cushions MS5 shown in Cherry. Four sizes of Seat Cushions and two Back Cushions for even more growth!
- 2 Headrest with Laterals MHL shown in Lime. See page 7.
- 79060159 Dark Gray, Small 79060359 Dark Gray, Large Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.
- **Standard Headrest** MHS shown in Cherry. See page 7.

Push Chair Sun Shade 79210000 Provides protection from the sun and a convenient window for monitoring.

## **Customize Your Size**



Small Shell (MS)

Large Shell (ML)

### Options for (MS) Shell

#### **Headrests**



7.5"W × 6.5"H



**Back Cushions** MS4



9"W × 11"H



#### **Seat Cushions**









### Options for (ML) Shell

#### Headrests



7.5"W × 6.5"H



### **Back Cushions**





#### **Seat Cushions**









### **Customize Your Size**

Fill in your custom sizes here for ease when ordering





**S**hells Choose a Shell size to determine your options









Seat

	Small MPS Seat Cushions				Large MPS Seat Cushions			
MPS Seat Specifications	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Wt.	80 lbs.				130 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16" 8" - 18"			- 18"				
Hip Flex Angle	-15° to +15°			-15° to +15°				
Seat Tilt Angle Forward	e Forward 10°			10°				
Seat Tilt Angle Back	k -30°		-30°					
Floor to Seat Ht.		14" -	28.5"		17" - 28.5"			

	Small MPS Back Cushions		Large MPS E	Back Cushions	
MPS Back Specifications	MS2 MS4		ML2	ML4	
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult	
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult	
Maximum User Wt.	80 lbs.		130	30 lbs.	
Back Width	8"	9"	10"	11"	
Back Height	8.5"	11"	11"	13"	
Top of Back Cushion	11" - 16" 15" - 20"		15" - 20"	19" - 22"	
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"	

MPS Tray	Small (79060155) Large (79060355)		
Dimensions	24"W × 17"D × 2"H		
Weight	8 lbs.		

MPS Hi-Low Seating System
79300105 Dark Gray, Small
79300305 Dark Gray, Large
Includes Seat Shell with Footrest,
8-point Harness, Hi-Low Base with
Casters (Must Choose Custom
Seating Cushions Separately)





		The second secon
Seat Cushions	Cherry	Lime
MS1	79050551	79050553
MS3	79050151	79050153
MS5	79050251	79050253
MS7	79050651	79050653
ML1	79050751	79050753
ML3	79050351	79050353
ML5	79050451	79050453
ML7	79050851	79050853

Back Cushions	Cherry	Lime
MS2	79120551	79120553
MS4	79120151	79120153
ML2*	79120751	79120753
ML4	79120351	79120353

Headrests	Cherry	Lime
MHS (Standard)	79100051	79100053
MHL (w/Laterals)	79110051	79110053

\*ML2 designed to fit both the Small and Large system, providing a wider Back Cushion option for Small size systems. Warranty: 5-Year Manufacturer's Warranty on cushions, and seat shell; 2-Year Manufacturer's Warranty on Hi-Low base, and strapping.



# **Soft-Touch® Sitters**

Supportive, Stylish Seating



# **Soft-Touch® Sitters**

Supportive, Stylish Seating

**Base Styles** Three styles to fit your needs



Sitter with Stationary Base



Sitter with Mobile Base



Floor Sitter Kit for Sizes 1,2,3

### **Features**



Sleek modern design provides up to 25° of tilt




Attachment straps for securing to a standard chair











# Soft-Touch<sup>®</sup> Sitters 5 year warranty





### **Applications** Sitters can be used in a variety of ways





Sitter Specifications	Size 1	Size 2	Size 3	Size 4	Size 5
Target Weight	20 - 30 lbs.	30 - 60 lbs.	40 - 80 lbs.	50 - 110 lbs.	90 - 200 lbs.
Target Height	30" - 35"	35" - 48"	40" - 56"	50" - 65"	54" - 74"
Target Age	1 - 2 years	2 - 6 years	6 - 9 years	9 - 14 years	13 years - Adult
Inside Seat Depth	7"	9.5"	11.5"	13.5"	16"
Inside Seat Width	8"	10"	12"	14"	16"
Seat Back Height	21"	26"	29"	35"	39"
Floor to Seat Ht. UPRIGHT on Stationary Base	10"	10"	9.75"	11.5"	12.5"
Floor to Seat Ht. TILTED on Stationary Base	10.5"	10.5"	10.5"	13"	14.25"
Floor to Seat Ht. UPRIGHT on Mobile Base	12.5"	13.25"	13"	14.75"	16"
Floor to Seat Ht. TILTED on Mobile Base	12.75"	13.75"	13.75"	16"	17.75"
Floor to Seat Ht. UPRIGHT on Floor Base	5"	5.5"	5.5"	N/A	N/A
Floor to Seat Ht. TILTED on Floor Base	4.5"	4.5"	5"	N/A	N/A
Overall Seat Depth	12"	14"	16"	19"	21"
Overall Seat Width	12"	14"	16"	18"	20"
Overall Seat Height	25"	29"	34"	40"	44"
Weight of Seat	6 lbs.	10 lbs.	13 lbs.	15 lbs.	19 lbs.

Base Specifications	Floor Base	Stationary (Small)	Stationary (Large)	Mobile Base (Small)	Mobile Base (Large)
Use with Sitter Sizes	1, 2, 3	1, 2, 3	4, 5	1, 2, 3	4, 5
Overall Width	19"	17.5"	19.25"	17.5"	19.25"
Overall Depth	22.5"	30"	35"	30"	35"
Overall Height	12.5"	16"	23"	19.5"	26.5"
Weight of Base	14 lbs.	12 lbs.	15 lbs.	14 lbs.	17 lbs.
Degree of Tilt	25°	25°	25°	25°	25°

Soft-Touch® Sitters (w/ Gray Harness)	Aqua	Dark Gray
Size 1 Sitter Only	77010104	77010109
Size 2 Sitter Only	77010204	77010209
Size 3 Sitter Only	77010304	77010309
Size 4 Sitter Only	77010404	77010409
Size 5 Sitter Only	77010504	77010509
Size 1 Sitter on Stationary Base	77110104	77110109
Size 2 Sitter on Stationary Base	77110204	77110209
Size 3 Sitter on Stationary Base	77110304	77110309
Size 4 Sitter on Stationary Base	77110404	77110409
Size 5 Sitter on Stationary Base	77110504	77110509
Size 1 Sitter on Mobile Base	77210104	77210109
Size 2 Sitter on Mobile Base	77210204	77210209
Size 3 Sitter on Mobile Base	77210304	77210309
Size 4 Sitter on Mobile Base	77210404	77210409
Size 5 Sitter on Mobile Base	77210504	77210509
Size 1 Floor Sitter	78010104	78010105
Size 2 Floor Sitter	78010204	78010209
Size 3 Floor Sitter	78010304	78010309

Optional Items	Part No.
Optional Floor Sitter Wedge (fits sizes 1-3 sitter)	78000009
Optional Small Stationary Base (fits sizes 1-3 sitter)	77100100
Optional Small Mobile Base (fits sizes 1-3 sitter)	77200100
Optional Headrest	77020000
Optional Large Stationary Base (fits sizes 4-5 sitter)	77100200
Optional Large Mobile Base (fits sizes 4-5 sitter)	77200200

# Mobile Activity Tray 2 year warranty



Portable and Versatile Tray for Everyday Use



- Two sizes: Short and Tall
- Height adjustable tray surface
- Angle adjustable tray surface
- Easy to move around a room
- Four tray surface color options







Mobile Activity Tray							
Mobile Activity Tray Short - Lilac	65000102	Mobile Activity Tray Tall - Lilac	65000302				
Mobile Activity Tray Short - Aqua	65000104	Mobile Activity Tray Tall - Aqua	65000304				
Mobile Activity Tray Short - Chocolate	65000105	Mobile Activity Tray Tall - Chocolate	65000305				
Mobile Activity Tray Short - Dark Gray	65000109	Mobile Activity Tray Tall - Dark Gray	65000309				

Mobile Activity Tray - Short Specifica	itions
Lowest Floor-to-tray surface	16" H
Highest Floor-to-tray surface	21" H
Weight	24lbs.
Tray Dimensions	28" W x 20" D

Mobile Activity Tray - Tall Specifications	
Lowest Floor-to-tray surface	20" H
Highest Floor-to-tray surface	28" H
Weight	29lbs.
Tray Dimensions	28" W x 20" D



### **Features**



- 1 Tray Tilt Angle 0 - 80 degrees.
- 4 Use with Wheelchair 28" wide base opening gives many wheelchair users access to the Mobile Activity Table.
- 2 Tray Height Adjustment Short: 16-21". Tall: 20-28".
- (5) Use with Soft-Touch® Sitters
  Designed to integrate with the
  Soft-Touch® Sitters for a variety
  of activities.
- 3 Colorful Options
  Vibrant color choices.
- 6 Use with Walker 28" wide base allows children with walkers access to the Mobile Activity Table from a standing position.



# Height Right<sup>™</sup> Chair 5 year warranty



### Simple Support that Grows

- Designed to grow with your child from 6 months up to 250 lbs
- Depth and height adjustable seat plate and foot plate
- NEW Early Pelvic Positioner offers contoured support starting from 6 months



3-Point Pelvic Belt Included.



Optional Soft-Touch® Seat and Back Cushions can be used for additional comfort.



Height Right™ Chair - 61000000	
Seat Plate	16.5"W × 9.25"D
Seat Plate Height Range (from floor)	16.5" - 24"
Foot Plate	11.5"W × 16.5"D
Foot Plate Height Range (from floor)	4" - 15"
Tray	19"W × 8"D
Tray Height to Floor	26"
Weight Capacity	250 lbs.
Overall	17"W × 21"D × 33"H

Early Pelvic Positioner Specifications	
Lateral Thigh Support Height	3.5"
Seat Width	11"
Seat Depth	9.5"

**Warning:** The Special Tomato Height Right™ Chair is as stable as a typical kitchen chair or baby highchair. However, if your child selfstimulates using vestibular movements (i.e. - rocking back and forth), this product might not be the best choice to meet your child's needs.





- 1 Soft-Touch® Seat Cushion 61040001 Cherry 61040002 Lilac 61040004 Aqua 61040009 Dark Gray Sold separately.
- Footcups 61060009 Dark Gray Only, Pair Sold separately.



**Warning:** The Early Pelvic Positioner is ALWAYS required when using the Tray for children UNDER the age of 3.

- 2 Soft-Touch® Back Cushion 61050001 Cherry 61050002 Lilac 61050004 Aqua 61050009 Dark Gray Sold separately.
- 5 Adjustable Foot Plate
  Shown without footcups

- 3 Early Pelvic Positioner 61070001 Cherry 61070002 Lilac 61070004 Aqua 61070009 Dark Gray Sold separately.
- 6 Wooden Tray with Cover 61010000 Sold separately.

Additional Tray Cover 61020000 Sold separately.



# Soft-Touch® Liners 5 year warranty



Upgrade Your Strollers and Chairs with Portable Positioning







Item	Cherry	Dark Gray	Inner Dimensions	Overall Dimensions
Size 1 Seat Cushion	72040101	72040109	8"W × 8"D × 3.5"H	10.5"D × 11"W
Size 2 Seat Cushion	72040201	72040209	10"W × 10"D × 4"H	13"D × 12"W
Size 3 Seat Cushion	72040301	72040309	12"W × 12"D × 4.5"H	15"D × 14"W
Size 1 Back Cushion*	72050101	72050109	7.5"W × 10"H	10"H × 10.5"W
Size 2 Back Cushion*	72050201	72050209	8.5"W × 13.5"H	13.5"H × 12.25"W
Size 3 Back Cushion*	72050301	72050309	10"W × 14.5"H	14.5"H × 13.25"W

<sup>\*</sup>Depth of Back Cushions are 1" at the narrowest point and 4" at the widest point



Please Note: Only Liner Sizes 1 and 2 can be inserted into the EIO Push Chair and Tomato Jogger.

# **5** year Soft-Touch® Booster Seats

Rise to Independence

- Secure 3-Point Pelvic Harness
- Provides height and pelvic stability
- Sacral and lateral thigh support





3-point harness.



Lightweight	and	portal	ole.
-------------	-----	--------	------

Soft <sup>-</sup> Touch® Booster Seat Specifications		
Lift of Seat off Chair	3.5"	
Seat Width	10"	
Seat Depth	9.5"	
Seat Back Height	4.5"	
Overall Height of Seat Cushion	7.75"	
Overall Dimensions	13.5"W × 13.5"D	



Please Note: The Special Tomato™ Soft-Touch® Booster has not been crash-tested and cannot be used as a car seat booster.




# 




**Convenient Comfort** 

**Elongated Potty Seat** 



Round Potty Seat





Item	Aqua
Round Seat Cushion	74010004
Elongated Seat Cushion	74020004





# **5**year Soft-Touch® Therapy Rolls & Wedges

Comfortable Therapeutic Support







Soft-Touch Rolls			
Diameter/ Length	Aqua Dark Gray		
3" × 8"	75000104	75000109	
4" × 24"	75000204	75000209	
6" × 24"	75000304	75000309	
8" × 24"	75000404	75000409	
8" × 36"	75000504	75000509	
10" × 36"	75000604	75000609	
12" × 36"	75000704	75000709	
12" × 48"	75000804	75000809	
14" × 48"	75000904	75000909	
16" × 48"	75001004	75001009	

Soft-Touch Wedges			
Elevation/ Width/Length	Incline	Aqua	Dark Gray
4" × 20" × 22"	10°	76000104	76000109
6" × 20" × 22"	16°	76000204	76000209
8" × 20" × 22"	20°	76000304	76000309
6" × 24" × 28"	12°	76000404	76000409
8" × 24" × 28"	16°	76000504	76000509
10" × 20" × 22"	26°	76000604	76000609
10" × 24" × 28"	20°	76000704	76000709
12" × 24" × 28"	24°	76000804	76000809

#### Rolls

- Great for postural therapy, balance, flexion and extension and neurodevelopment
- Solid Core provides firm platform (except 75000104, 75000205, 75000204 and 75000205)
- Stand-On-End™ Design for security and stability when not in use

#### Wedges

- Ideal for adaptive positioning, balance, and muscle tone control
- Also recommended for dynamic movement patterns and fine and gross motor activities
- Stand-On-End™ Design for security and stability when not in use
- Integrated Soft-Touch® handle for easy carrying



# EIO Push Chair 2 year warranty



### Affordable Chair and Positioning

- Designed for children from early intervention to pre-adolescent ages
- Easy to fold and transport
- Lightweight aluminum frame
- Calf panel and seat back lay your child flat
- Standard features: see next page

Strallar Only Soat to Foot Wall
Stroller Only Seat to Foot Well
Liners Size 1 and 2 Seat to Foot Well
Sitter Size 1 Seat to Foot Well
itter Size 2 Seat to Foot Well
Open Chair Dimensions
Folded Chair Dimensions
Handle Height to Floor
Seat
Back Height
Max. User Weight
Stroller Weight
Includes 2 foot-block





- 1 Size 1 and 2 Soft-Touch® Liners
  Attach to the stroller, with color choice
  of Cherry, Lilac, Aqua or Dark Gray.
  Sold separately, see page 16.
- 4 Shopping Basket
  10 lb. weight capacity.
- Sun Shade
  Adjustable push handle, padded headrest and removable sun canopy with window, which includes a large pocket.

- Size 1 and 2 Soft-Touch® Sitters Attach to the stroller, with color choice of Lilac, Aqua or Dark Gray. Sold separately, see page 9.
- 5 Angle Adjustable Seat Back (25° to 75°) and foot well (+15° to -55°).
- Rain Cover 63030000 Sold separately.



- 3 Adjustable 8-point Harness
  5-Point comfort trunk harness with padded shoulders and chest strap.
  Also includes a padded 3-point pelvic harness for hip and lateral support.
- Solid Filled Wheels Solid filled 9" front and 11" rear wheels require no air. Swivel front wheels with lock-out feature add security on rough terrain. Includes one-touch parking brake.

# Jogger 2 year warranty



### Supportive Style On-the-Go

- Designed for parents and families who want to take long strolls through the city or through rough trails
- Easy to fold and transport
- Lightweight aluminum frame
- Includes 8-point trunk / pelvic harness system







- 1 Size 1 and 2 Soft-Touch® Liners
  Attach to the stroller, with color choice
  of Cherry, Lilac, Aqua or Dark Gray.
  Sold separately, see page 16.
- Pneumatic Wheels Front 12" pneumatic swivel wheel with lock-out and quick-release 12" rear pneumatic wheels.
- Medical Necessities Bag 62010000 Sold separately.

- Size 1 and 2 Soft-Touch® Sitters Attach to the stroller, with color choice of Lilac, Aqua or Dark Gray. Sold separately, see page 9.
- **One-step, compact folding design** Canopy folds with Jogger.
- Bead Support 62020000 Sold separately.
- Plomoto moorings mediquip

Added Comfort
 Padded upholstery with

Padded upholstery with reclining back support, functional rear suspension and foot well with side guard.

- Sun Shade Removable canopy with window and height adjustable handle.
- **Rain Cover**62030000Sold separately.

# **The Special Tomato Story**



My dad and mom had four children. The first three were healthy and regular people. Then I came along; a baby with a difference—Down syndrome.

One day it was bothering me that I was having a hard time doing my school work. My dad shared a story about his tomato garden. Once, at harvest time, a differently shaped tomato got his attention. It had a large, bright red bump that made my dad want to pick it up and study it because it was so unique and looked appetizing. At dinner, that tomato was just as juicy as the others.

Even though I'm different—like the tomato—as a person with Down syndrome, my lovable and outgoing personality and strong character make me just as special as everyone else. I guess I am the special "tomato" in our family's garden!









**Customer Service Freephone 0800 031 6571** 

### **Moorings Mediquip UK Ltd**

Head Office and Accounts
51 Slaght Road,
Ballymena,
Co. Antrim
BT24 2JH

Tel: 028 2563 2777 Fax: 028 2563 2772 Email: sales@mooringsmediquip.com www.mooringsmediquip.com

www.specialtomatouk.co.uk