

## INTRODUCING...

standard 7176-19.

Ibis the complete tilt-in-space modular seating wheelchair, the base is combined with the award winning Sedeo Pro seating system enabling the Ibis to be one of the most adjustable wheelchairs in its class!

Introducing the tilt-in-space chair which can be simply adjusted through the use of one tool, with a wide range of adjustment available in seat width, depth and height. With 30° tilt the lbis has been designed specifically to meet the requirements of a wide range of users' needs, from those requiring comfort seating through to conditions such as MND, CVA and Muscular Dystrophy.

# **EXTREME ADJUSTMENT...**



## TAKE A SEAT...

Utilising the award winning Sedeo seating system the Ibis enables the chair to be superior in adjustability, comfort and support.

The lbis comes with an adjustable flexible back support ensuring the optimal postural management position is always achieved as the chair can be configured exactly to the user's specific requirements, in either a fixed or dynamic set up.

The thoracic supports of the 14cm & 21cm deep contour active back can be easily positioned in either a symmetrical position to promote mid-line sitting, or asymmetrically to support preferred angles or Scoliosis.

Kyphotic and Lordotic spinal presentations can be supported using the tension adjustable straps contained within the back.

The simple yet effective system ensures that the user is positioned optimally for their level of functional ability, whilst at the same time helping to prevent further physical deterioration.





# YOUR CHAIR YOUR CHOICE...

Ibis Pro is supremely set to meet the broadest range of requirements.

Careful selection of accessories and fittings have created a chair that, with the choice of back support and cushion will meet the user's size and posture requirements. The resulting seating system will then provide the user with a seating position that allows for them to maximise their functionality with minimum effort whilst reducing potential damage to themselves.

If a standard configuration is insufficient there are a broad range of accessories to meet those greater needs; negative tilt to aid with standing transfers or to promote active extension, thigh supports to help control abduction, manually elevating legrests for those users with differing flexibility or even a Power Support unit to make life easier and less stressful for both the user and carer.

The lbis can easily accommodate third party seating systems which ensures that the lbis is a perfect choice for many users.



For further information on the full specification, options and accessories, please refer to the order form on our website.

All information is subject to change without notification. Please consult Sunrise Medical with any gueries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively it is available on request in large text.

















#### Sunrise Medical Limited

Thorns Road, Brierley Hill, West Midlands, DY5 2LD

Tel: +44 (0) 845 605 66 88

E-mail: enquiries@sunmed.co.uk

www.SunriseMedical.co.uk