

E|F55

Retail Order Form

**UK
CA
CE**


ORDER FORM

Billing Address Or Customer-ID		Shipping Address		<input type="checkbox"/> Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
<input type="checkbox"/> Order <input type="checkbox"/> Calculation <input type="checkbox"/> Built-4-Me®		Max. user weight: 110 kg		
Purchase Order No.: <input type="text"/>		User Weight: <input type="text"/> kg		
Art. No.: 47600000				
Order Date: Day <input type="text"/> Month <input type="text"/> Year <input type="text"/>				

Technical Data Overview

Max. user weight	110 kg	Max. speed forward	6 kph; optional 10 kph, 15 kph or 20kph
Front frame width	260 – 480 mm ⁽¹⁾	Wheel size	8,5" air-chambered solid-tyre / 14" pneumatic tyre
Battery weight	2 kg	Battery type	Li-Ion, 36V, 8.1 Ah
Total weight w/o battery	11 kg	Brake type	Twin-brake system: mechanical disc-brake and separate electric brake
Total weight with battery	13 kg	Operating Voltage	36 Volt
Colour	Black	Range of action	25 km max.
Colour accents (included in delivery)	Black, red, blue, orange	Climbing power	up to 11% (6°)

Base		£
<input type="checkbox"/>	Pull-device F55 400 W (max.) motor, 8.1 Ah Li-ion battery, forward- and reverse gear, cruise-control, dual-brake-system (disc-brake, e-brake), mudguard, front-light and colour-kit included.	0
Drive wheels		£
<input type="checkbox"/>	EMF070085 Drive wheel, 8,5", solid tyre	0
<input type="checkbox"/>	EMF070086 Drive wheel, 14", pneumatic tyre	0
Ordering/Bundling		£
<input type="checkbox"/>	EMF090213 Without wheelchair order	2,995
<input type="checkbox"/>	EMF090217 With wheelchair order, not mounted to the wheelchair (1)	2,995
Docking		£
<input type="checkbox"/>	EMF010027 Docking for fixed-front wheelchairs	0
<input type="checkbox"/>	EMF010028 Docking for swing-away wheelchairs	0
Docking-Clamps		£
<input type="checkbox"/>	EMF010042 Docking-clamps for Ki Mobility Rogue (3) (4)	0
<input type="checkbox"/>	EMF010041 Docking-clamps for oval tubes (Progeo Joker / Ego / Tekna) (3) (4)	0
<input type="checkbox"/>	EMF010029 Docking-clamps for round tube frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm & 30 mm tubes) (3) (4)	0
<input type="checkbox"/>	EMF010030 Docking-clamps for oval tubes (Quickie Xenon ² FF, Helium, Nitrum) (3)	0
<input type="checkbox"/>	EMF010031 Docking-clamps for Quickie Xenon ² SA (3)	0
Adjustability		£
<input type="checkbox"/>	Complete easy-adapt kit included to fit F55 to front frame widths from 260 – 480 mm (5). The F55 is fully adjustable to all seat heights and frame angles.	0
Speed Setting		£
<input type="checkbox"/>	EMF000025 6 kph	0
<input type="checkbox"/>	EMF000027 10 kph (6)	0
<input type="checkbox"/>	EMF000028 15 kph (6)	0
<input type="checkbox"/>	EMF000167 20kph (6)	0
Display		£
<input type="checkbox"/>	Display: Battery status, Odometer, trip counter, avg. speed, real time power consumption	0
Accessories		£
<input checked="" type="checkbox"/>	LED front light	0
<input checked="" type="checkbox"/>	Mudguard for front wheel	0
<input checked="" type="checkbox"/>	Colour-kit included: mudguard and sidebars in 4 colours	0
<input type="checkbox"/>	EMF010032 Kickstand: mounts to the sub-frame	185
<input checked="" type="checkbox"/>	EMF010033 Padding for sub-frame (enclosed)	0
<input type="checkbox"/>	EMF160010 Additional Li-Ion battery, Standard: 36V, 8.1 Ah, 291 Wh (approx. 25 km range)	500
<input type="checkbox"/>	EMF110155 Additional charger	130

(1) Wheelchair order on separate order form. Separate shipment of wheelchair and F55.

(2) Wheelchair-order on separate order form. Shipment of readily mounted wheelchair/F55-unit

(3) Always check installation space for the docking clamps before ordering. Depending on configuration and option, a particular configuration of a chair might not offer enough space.

(4) Always check combination agreements for non-Sunrise wheelchairs.


(5) Frame width measured from center frame tube to center frame tube at the mounting position. Use only the measured dimension and not the seat width to check compatibility.

(6) Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via e-mail with a signed waiver (next page). Max. speed-settings can only be altered by Sunrise Medical.

Accessories



LED-light, Fender (accent colour red) and disc-brake (pictured on EMF070085 Drive wheel 8,5", solid tyre)



// EMF010032 Kickstand: mounted to the sub-frame



Colour-kit (mudguard) included in delivery



Colour-kit (sidebars) included in delivery (accent colour blue)

Letter of Acknowledgement

We herewith advise, that using the Empulse F55 with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the Empulse F55 with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Empulse F55 speed settings higher than 6 kph (10/15/20 kph). When you're travelling with your F55, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your F55 on private property.

I herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.

Date

Customer Signature

Day

Month

Year

Comments

F55

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day

Month

Year

Send E-Mail