

Quick To Connect:

The patented One Touch Clamps enable easy, fast and effortless docking to the wheelchair. When undocking, nothing is left on the chair.



Powerful To Drive:

High performance motor and a high battery capacity will deliver up to 50 km of range, for worry-free driving fun. Available with optional speed settings of 6, 10, 15, and now 20kph.



Easy To Use:

The StreetJet was designed to be as easy as possible to operate, meaning even users with limited hand function can experience the benefits.

Book a demo today.



Technical Data

Range of action	up to 50 km ¹
Max. user weight	120 kg
Max. speed forward	6 kph; optional 10, 15 or 20 kph
Motor power (max power)	12" 250 W (630 W) / 16" 250 W (750 W)
Max safe slope	6°
Front frame width	220 - 425 mm ²
Battery weight	2.6 kg
Total weight w/o battery	from 12,9 kg ³
Total weight with battery	from 15.5 kg ³
Operating voltage	36 Volt
Battery type	Lithium-ion battery, 36 V 14 Ah, intelligent battery management system
Charger	36V, 2A; optional 4A quick
Wheel size	12" pneumatic tyre / 16" pneumatic tyre
Brake type	Twin-brake system: mechanical disc brake and separate proportional electric brake. Optional Triple-brake system: dual mechanical brake and separate electric brake

Key Features

Patented, Fast And Effortless Docking:

The patented One Touch Clamps were developed to enable easy, fast and effortless docking to your wheelchair. Designed to fit a wide range of frame shapes, the StreetJet doesn't require any additional accessories on your wheelchair, which would interfere and add weight when not in use.

More Power For Extended Tours:

Thanks to StreetJet's power you can satisfy your longing for adventure and experience a new sense of freedom and independence. Take advantage of the extra power the motor provides to go longer distances and don't let tired shoulders dictate where you go.

New 20kph Max Speed Option:

Optional speed settings of 6, 10, 15, and now 20kph make the StreetJet ideal for any adventure. You can also set cruise control, allowing you to maintain speed for longer periods without using the throttle.

The Right Wheel Size For You:

Looking to cruise around town? The 12" wheel offers a shorter base length and smaller turning circle - perfect for urban environments and for transportation. Alternatively, if you prefer more outdoor adventures, the 16" wheel offers increased comfort for tackling difficult terrain - made for the great outdoors.

Customise Your StreetJet

Set the tone of your StreetJet with 20 frame and fork colours.



Certified according to all EN12184 applicable standards

(1) According to ISO 7176 - 4, depending on the ground surface, the weight of the user, driving behaviour, temperature and terrain. (2) Measured inside tube to inside tube. Use only the measured dimension and not the seat width to check compatibility. (3) Depending on product configuration

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines.
Manufacturer: Zhejjang Innuovo Rehabilitation, Devices Co. Ltd, No. 1896 Industry Road, Hengdian Movie Zone, Dongyang, Zhejjang, China
In line with ongoing product improvement, Innuovo reserves the right to make changes at any time, and without notice. However, any changes to information provided for existing users shall be clearly communicated if they are safety critical. Not all features and options are compatible with all configurations. All dimensions are approximate and may be subject to change.

