



With the choice of a RIGID frame for performance or a FOLDING frame for easier transfers, it is easy to find a Life to suit you.

Ouickie Life, folding fixed front



Ouickie Life, folding swing-away



Quickie Life R, rigid fixed front



Quickie Life R, rigid swing-away







Rigid, durable and robust the Life R is designed to enhance your life now and for years into the future.

Available as fixed front or swing-away, the new Life R has a modular design which means the front of the frame can be changed to increase seat depth, or alter the frame type. This means it's perfect for growing teenagers and makes it easy to refurbish.



Whatever the day throws at you, the energy efficient design of the Life R means you can focus on what's important to you.







Easy Life

In a business meeting this morning? Preparing for the great outdoors this afternoon? Whatever your plans your Quickie Life is ready.

If transportation and a small folding package are key, choose the Life folding. With full adjustability and configurability the Quickie Life is perfect for those who lead an active lifestyle.

But why stop there? Personalise your Life by selecting from a wide range of side guards, armrests, hangers, footrests, wheels and handrims. Or define your style with a choice of 10 frame colours and anodized parts.



Combine the modern design with new frame colours and anodized parts for a stylish chair to match your life.





Configure your Life with a choice of frame styles, angles and lengths.





Seat height adjustable in seven positions

With a camber from 0° to 4° and four centre of gravity positions (from the very active 8.3cm to the more steady 2.3cm).



Active frame (optional)

Compact frame (85°) with 2cm inset for more active users.



Simple and quick to adjust fork angle

With increased angle range to ensure optimal driving performance.



Hand bike axle

Increased stability to use the chair in combination with a hand bike.



Angle adjustable back (optional)

Lightweight and simple to adjust. Provides 20° of angle adjustment from -12° to +8° in 4° increments.



Lightweight platform footrest

Angle and depth adjustable.

Available in composite and aluminium.



Lightweight height adjustable single post armrest

Provides easy transfer and good arm support with a low transport weight.



Active swing-away hangers (optional on swing-away frame)

Offers a modern and active look. Available in 70° and 80°.



Whatever your lifestyle, Quickie Life folding can be adjusted and configured to suit. Choose from a fixed front lightweight frame for a more active look, or the swing-away version for even easier folding and transportation, both now available with new lightweight options.





NEW Aluminium lightweight hangers in 70° and 80°. 13% lighter than Quickie standard hangers.



NEW Lightweight height adjustable single post armrest
Provides easy transfer and good arm support with a low

transport weight.



NEW Lightweight anti-tip for optimal driving performance with 480g weight reduction.



Angle adjustable back (optional)
Lightweight and simple to adjust. Provides angle adjustment from -12° to +12° in 4° increments.



Active frame
Compact frame (85°) with a 2cm inset for more active users. Available with or without a 2cm inset.



NEW folding fixationOperated with one hand, keeps the chair securely folded.



Flip-up platform footrest Lightweight, angle and depth adjustable.



Seat height adjustable in seven positions
With a camber from 0° to 4°, and 4 centre of gravity
positions (from the very active 8.3 cm to the more stable
2.3 cm).



Technical data sheet



| | Life fixed front | Life swing-away | Life R fixed front | Life R swing-away |
|--------------------------------------|-------------------------------|--------------------|------------------------------|--------------------|
| Overall width | Seat width +220mm | SW + 220mm | SW + 220mm | SW + 220mm |
| Seat width | 360 – 500mm | 360 – 500mm | 340 – 480mm | 340 — 480mm |
| Back height | 250 – 475mm | 250 – 475mm | 250 – 450mm | 250 – 450mm |
| Seat depth | 360 – 500mm | 360 – 500mm | 360 – 500mm | 360 — 500mm |
| Rear seat height | 400 – 500mm | 400 – 500mm | 360 – 450mm | 360 – 450mm |
| Front seat height | 400 – 530mm | 400 – 530mm | 390 – 500mm | 390 – 500mm |
| Back angle | -12° to +12° | -12° to +12° | -12° to +8° | -12° to +8° |
| Frame front end angle | 75°, 85°, 85° with 20mm inset | 80° and 70° | 75°, 85°, 85° with 2cm inset | 80° and 70° |
| Individual frame colours | 32 | 32 | 32 | 32 |
| Individual sling colours | 1 | 1 | | 1000 |
| Fork and axle plate anodized colours | 3 | 3 | 3 | 3 |
| Total chair weight, from | 11.5kg | 12kg | 10.7kg | 11.3kg |
| Maximum user weight | 125kg (19.6 stone) | 125kg (19.6 stone) | 125kg (19.6 stone) | 125kg (19.6 stone) |



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Life family and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairs



Sunrise Medical, Thorns Road, Brierley Hill, West Midlands, DY5 2LD Tel: +44 (0) 845 605 66 88 Fax: +44 (0) 845 605 66 89 www.SunriseMedical.co.uk