



SALSA Mini²

RED LINE

MINI SIZE

MAXIMUM PERFORMANCE



NEW SALSA M² MINI *RED LINE*

Introducing the Salsa M² Mini *Red Line*
- when maximum performance and size
matter.

Pure *Speed* - Red line the speed gauge with
high-torque 10kph motors and traction control technology
ensuring you're always in control.

Pure *Agility* - Access all areas with a super slim 52cm
base, it's the narrowest mid-wheel drive performance
powered wheelchair in the world!

Pure *Style* - Beautifully designed, a sleek matt black
chassis and base are perfectly complemented by red trim
on the upholstery, springs and details.

Pure *Luxury* - Take a seat and experience true luxury
with our top of the range seating, the JAY Comfort.

DON'T COMPROMISE.



PURE **INDOOR** AGILITY

With its slim base of just 52cm, the Salsa M² Mini *Red Line* is the narrowest high-performance mid-wheel-drive powered wheelchair in the world, providing you with unparalleled indoor agility.

Restaurants, shops, or at home - condensed layouts are easy to negotiate with the super-narrow Red Line.

Access through floor lifts and hoists - where other powered wheelchairs might be too big or heavy. The Red line is narrower and lighter (from just 100kg) than any other mid-wheel-drive performance powerchair on the market.

Be in the driving seat - Easily access Wheelchair Adaptive Vehicles (WAVs) and comfortably sit under the steering wheel with the Salsa M² Mini Red Line's ultra-low seat to floor height, anti-pitch technology and narrow dimensions.





10kph motors



LED lights and indicators




6 wheel independent suspension

PURE OUTDOOR PERFORMANCE


Highly-advanced, proven technology that's finely engineered for outdoor performance.

10kph motors make easy work of slopes of up to 8° and kerbs of up to 3", whilst integrated Gyroscopic technology works like traction control keeping you on the straight and narrow by auto-correcting any deviations on gradients and Quickie's patented all-wheel suspension absorbs all those bumps to give you a smooth and stable ride over slopes and challenging terrain.

With the same mid-wheel-drive technology as the original Salsa M² Mini your weight is positioned directly over the drive wheel – maximising traction and providing a safe driving position.



Without Gyro Control Module



With Gyro Control Module



The Salsa M² Mini Red Line is transit tested in accordance with ISO 7176-19 and approved for transportation in a vehicle. The tests are conducted with a four-point tie-down system.





PURE STYLE

Beautiful, sleek and refined design means the *Red Line* has the most minimalistic and modern style.

Experience a new level of style with the Quickie Salsa M² Mini Red Line. The matt black base, chassis, castor arms and upholstery are perfectly complemented by the red trimmed upholstery, springs and detailing.

PURE LUXURY

Take a seat on the NEW Quickie Salsa M² Mini *Red Line* and experience true luxury.

With our top of the range seating solution the JAY Comfort, enjoy a supportive, comfortable and speedy ride tailored to you.



PURE CONNECTIVITY

Salsa M² Mini *Red Line* boasts the option of the NEW LCD Advanced Joystick. Offering the ability to connect with your TV through Infra-Red or up to two bluetooth devices and one iOS device.

// BLUETOOTH & I-DEVICES

Connect to a range of technology, including PCs, laptops, phones and tablets with your Joystick, exclusively available on the Salsa M² Mini Red Line. It allows you to connect and take control of two Bluetooth and one iOS device.

// INFRA-RED

Control your home environment with the integrated Infra-Red technology, allowing you to connect to your TV, DVD and other devices giving you total control from the comfort of your chair.

// QUICKIE PATENTED ASSIGNABLE BUTTONS

Get super-quick access to commonly used commands by assigning them to any button on your joystick and assign second (long press) functions for total control at your fingertips.



TECHNICAL DATA

SALSA M ² MINI RED LINE			
Max. user weight:	140kg (22 Stone)	Battery size:	36Ah
Front & rear castor wheels:	180mm (7") Solid	Max. dynamic stability:	8° (Depending on options)
Drive wheels:	320mm (12.5")	Seat width & overall width: adult seat with standard armpads	510mm (20") - 630mm (25") 480mm (19") - 605mm (24") 455mm (18") - 580mm (23") 430mm (17") - 555mm (22") 405mm (16") - 530mm (21")
Seat width:	410 - 510mm (16" - 20")	Seat Depth:	410 - 510mm (16" - 20")
Min turning radius:	700mm (27.6") - with uni tyre hangers	Max. kerb climbing:	70mm (3")
Max. speed:	10kph (6mph)	Colour choice:	Matt Black
Max. range (ISO 7176-4):	23km (14 miles)	Transport (ISO 7176-19):	4 point tie-down & Dahl docking system
Battery box to ground clearance:	65mm (2.5")	Chair weight:	97kg - 128kg (Depending on options)



For further information on the full specification, options and accessories, please refer to the order form on our website. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively it is available on request in large text.

Sunrise Medical Limited
 Thorns Road, Brierley Hill, West Midlands, DY5 2LD
 Tel: +44 (0) 845 605 66 88
 E-mail: enquiries@sunmed.co.uk
www.SunriseMedical.co.uk