2017/18

ottobock.

Foot and leg
positioning
components



		Order no.:
Provide estimate only	□Order	
Billing address/customer number		
Delivery address		
	1	
Sales person	Name	
	 Date Sig	gnature
By placing the order, you agree to comply with accessories and products of other manufacture product, in addition to the instructions for use products. Also you confirm to have read, under	ers when combining /service manuals of t	them with an Ottobock he respective Ottobock

accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Foot straps and fo	ot caps					
Article no.		Figure	Designation	Dimensions in mm	Quantity	
476P52=SK030			Foot strap*, size XS (pair)	155–200		£157.00
476P52=SK031			Foot strap*, size S (pair)	160–215		£157.00
476P52=SK032			Foot strap*, size M (pair)	175–230		£157.00
476P52=SK033			Foot strap*, size L (pair)	240–295		£157.00
476P52=SK034			Foot strap*, size XL (pair)	295–340		£157.00
476B52=SK010		I S	Foot caps* – "small"	L: 190, W: 89		£87.00
476B52=SK015			Foot caps* – "medium"	L: 230, W: 114		£87.00
476B52=SK020			Foot caps* – "large"	L: 255, W: 140		£87.00
* For appropriate for	ot positic	ning on wheelchair	foot supports/ foot support plates, als	so can be used with Reha-Buggies	or mobility bases.	

Otto Bock Mobility Solutions Plc - 32 Parsonage Road, Englefield Green, Egham, Surrey, TW20 0LD - Tel 01784 744900 - Fax 01784 744901 - Email bockuk@ottobock.com, Web www.ottobock.co.uk

Technical modifications and printing errors reserved

Foot buckets and	calf sup	ports				
Article no.		Figure	Designation	Dimensions in mm	Quantity	
440R9=1-7			Foot bucket (large)	H: 325, W: 185, L: 255		£74.00
440R9=2-7		a de la compañía de	Foot bucket (medium)	H: 220, W: 150, L: 255		£74.00
440C2=1-7			Foot bucket padding (large)	H: 325, W: 185, L: 255		£38.00
440C2=2-7		S. D.	Foot bucket padding (small)	H: 220, W: 150, L: 255		£38.00
440P1=1-7			Calf support	B: ~167, H: ~230		£42.00
440P2=1-7*			Calf support	B: ~164, H: ~230		£33.00
* Note: Possible in c	ombinati	on with S.T.S. moun	ting kit (434F3*) and S.T.S. clamps (434H1*)). See the following table.		
S.T.S. mounting ki	its and cl	lamps				
Article no.		Figure	Designation	Dimensions in mm	Quantity	
434F3=L			S.T.S. mounting kit, left, offset,	A: 110, B: 95, C: 45		£32.00

434F3=L			S.T.S. mounting kit, left, offset, chrome	A: 110, B: 95, C: 45	£32.00
434F3=R			S.T.S. mounting kit, right, offset, chrome	A: 110, B: 95, C: 45	£32.00
434F3=L-7			S.T.S. mounting kit, left, offset, black	A: 110, B: 95, C: 45	£32.00
434F3=R-7			S.T.S. mounting kit, right, offset, black	A: 110, B: 95, C: 45	£32.00
434H1=3/4			S.T.S. clamp, chrome	Tube Ø 19.0	£13.00
434H1=3/4-7			S.T.S. clamp, black	Tube Ø 19.0	£13.00
434H1=7/8			S.T.S. clamp, chrome	Tube Ø 22.2	£13.00
434H1=7/8-7			S.T.S. clamp, black	Tube Ø 22.2	£13.00
434H1=1			S.T.S. clamp, chrome	Tube Ø 25	£13.00
434H1=1-7			S.T.S. clamp, black	Tube Ø 25	£13.00

Once printed this document is an uncontrolled version. Always check the download centre on our website ensure you have the most up to date version of this document.