

Invacare® | Softform® Premier HS

A clinically tested mattress with integral heel protection



Integral heel protection

A clinically tested, high specification foam mattress suitable for 'Very High Risk' patients.

The *Invacare Softform Premier HS* is a high specification foam pressure redistributing mattress designed to meet the demands of Community and Acute healthcare environments. Practical and durable, this static mattress delivers exceptional levels of comfort and pressure redistribution. Designed with a heel slope in the foot end, the mattress offers additional protection to the vulnerable heel area.

The mattress features Invacare's 'Strikethrough Resistant Technology'2 cover. For more information, please see brochure number 1505337.

The Invacare Softform Premier HS features a unique

contoured foam insert, constructed from 'teardrop

cut' castellated cells, which conform independently

and friction forces. In addition it provides excellent

redistribution. The mattress insert has a non-glued

construction, which allows for excellent conformity

to body movement and position to reduce shear

comfort capabilities and improved pressure

when profiled on electric bed frames.

The comfort of 'high specification' foam

Additional protection and redistribution for heels

Designed to address an important risk area in pressure ulcer development, the **Softform Premier HS** incorporates a dedicated heel protection area. The lower third of the mattress has smaller castellations to offer pressure redistribution to the bony prominence of the heel. In addition, the foam has a chamfered design which slopes toward the foot end to effectively support the lower leg whilst offering additional protection for the heel. The tear drop design of the castellations also offers protection from shear forces when the bed is profiled. 3

The Invacare Softform Premier HS is based on the popular Softform Premier mattress. For further information on the features and benefits of the *Invacare Softform Premier* mattress, please refer to brochure number 1505337.

Technical data

For more comprehensive pre-sales information about the Softform Premier, including the product's user manual, please visit www.invacare.co.uk.



Length 660 - 900 mm 1410 - 2100 mm



152 mm





39 Stone



Manufactured to European Standards: BS EN 597-1 Cig BS EN 597-2 Match. Foam Interior: BS 5852: Crib Source 5; Cover BS7175: Crib Source 5.

Softform Premier HS

Violet

Colours

Warranty

10 years foam 5 years cover



Single sizes only

Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement), Informal assessment is considered to be of greater importance and clinical value.

When profiling a patient in bed, it is considered best practice to elevate the knee break prior to elevating the backrest section or if available to use the auto profile function, this helps to reduce the risk of the patient migrating downward in bed and may reduce the risk of shear forces.

Invacare Limited Pencoed Technology Park - Pencoed Bridgend CF35 5AQ - United Kingdom Tel: +44 1656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk

Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS -County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie







